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PARAMETERS							
N°	@	Test	ACUR	Method	End-Point	Direction	Ascend
Unit	mg/dL	Decimal	1	Prim. Wavelength	510	Sec. Wavelength	700
Sample Volume	5	R1 Volume	200	R2 Volume	0	Line Limit	N/A
				Response			
Incubation	--	Check Antigen	No	Lower	--		
Reaction	0 -- 25	Substrate	N/A	Upper	--		
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0.0	Print Number	@
Upper	@	Upper	@	Upper	25.0		@
Calibration							
Method	Linear	Std. 0	0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	ALB	Method	End-Point	Direction	Ascend
Unit	g/dL	Decimal	1	Prim. Wavelength	630	Sec. Wavelength	0
Sample Volume	3	R1 Volume	300	R2 Volume	0	Line Limit	N/A
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	0 -- 10	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0.1	Print Number	@
Upper	@	Upper	@	Upper	6.0		@
Calibration							
Method	Linear	Std. 0	0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	AEO	Method	Fixed-Time	Direction	Ascend
Unit	UI/mL	Decimal	0	Prim. Wavelength	546	Sec. Wavelength	0
Sample Volume	3	R1 Volume	300	R2 Volume	0	Line Limit	N/A
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	3 -- 13	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	1	Print Number	@
Upper	@	Upper	@	Upper	800		
Calibration							
Method	Linear	Std. 0	0.0 (NaCl)				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	AEO	Method	Fixed-Time	Direction	Ascend
Unit	UI/mL	Decimal	0	Prim. Wavelength	546	Sec. Wavelength	0
Sample Volume	3	R1 Volume	240	R2 Volume	60	Line Limit	N/A
				Response			
Incubation	3	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	3 -- 14	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	1	Print Number	@
Upper	@	Upper	@	Upper	800		
Calibration							
Method	Linear	Std. 0	0.0 (NaCl)				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	ALT	Method	Kinetic	Direction	Descend
Unit	U/L	Decimal	0	Prim. Wavelength	340	Sec. Wavelength	670
Sample Volume	20	R1 Volume	200	R2 Volume	0	Line Limit	20
Incubation	--	Check Antigen	No	Response			
Reaction	3 -- 10	Substrate	N/A	Lower	--	Sample Volume	@
				Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0	Print Number	@
Upper	@	Upper	@	Upper	400		@
Calibration							
Method	Linear	Std. 0	0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	AMIL	Method	Kinetic	Direction	Ascend
Unit	U/L	Decimal	0	Prim. Wavelength	405	Sec. Wavelength	670
Sample Volume	4	R1 Volume	200	R2 Volume	0	Line Limit	20
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	5 -- 20	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0	Print Number	@
Upper	@	Upper	@	Upper	2000		@
Calibration							
Method	Linear	Std. 0	0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	AST	Method	Kinetic	Direction	Descend
Unit	U/L	Decimal	0	Prim. Wavelength	340	Sec. Wavelength	670
Sample Volume	20	R1 Volume	200	R2 Volume	0	Line Limit	20
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	3 -- 10	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0	Print Number	@
Upper	@	Upper	@	Upper	400		@
Calibration							
Method	Linear	Std. 0	0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	Ca Ars	Method	End Point	Direction	Ascend
Unit	mg/dL	Decimal	1	Prim. Wavelength	670	Sec. Wavelength	0
Sample Volume	2.5	R1 Volume	250	R2 Volume	0	Line Limit	N/A
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	0 -- 7	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0	Print Number	@
Upper	@	Upper	@	Upper	17		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	Ca CRE	Method	End Point	Direction	Ascend
Unit	mg/dL	Decimal	1	Prim. Wavelength	578	Sec. Wavelength	670
Sample Volume	4	R1 Volume	200	R2 Volume	0	Line Limit	N/A
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	0 -- 7	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0	Print Number	@
Upper	@	Upper	@	Upper	16		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	CK	Method	Kinetic	Direction	Ascend
Unit	U/L	Decimal	0	Prim. Wavelength	340	Sec. Wavelength	670
Sample Volume	4	R1 Volume	200	R2 Volume	0	Line Limit	20
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	10 -- 20	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0	Print Number	@
Upper	@	Upper	@	Upper	2000		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	CK MB	Method	Fixed Time	Direction	Ascend
Unit	U/L	Decimal	0	Prim. Wavelength	340	Sec. Wavelength	670
Sample Volume	10	R1 Volume	200	R2 Volume	0	Line Limit	20
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	25 -- 50	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0	Print Number	@
Upper	@	Upper	@	Upper	500		
Calibration							
Method	linear						

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PARAMETERS							
N°	@	Test	COL	Method	End Point	Direction	Ascend
Unit	mg/dL	Decimal	0	Prim. Wavelength	510	Sec. Wavelength	700
Sample Volume	3	R1 Volume	300	R2 Volume	0	Line Limit	N/A
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	0 -- 50	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	1	Print Number	@
Upper	@	Upper	@	Upper	500		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	HDL	Method	End Point	Direction	Ascend
Unit	mg/dL	Decimal	0	Prim. Wavelength	510	Sec. Wavelength	0
Sample Volume	20	R1 Volume	200	R2 Volume	0	Line Limit	N/A
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	0 -- 50	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	2	Print Number	@
Upper	@	Upper	@	Upper	150		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	40				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	CHE	Method	Kinetic	Direction	Descend
Unit	U/L	Decimal	0	Prim. Wavelength	405	Sec. Wavelength	0
Sample Volume	4	R1 Volume	200	R2 Volume	50	Line Limit	20
				Response			
Incubation	15	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	10 -- 20	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	123	Print Number	@
Upper	@	Upper	@	Upper	20000		
Calibration							
Method							
Factor							

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PARAMETERS							
N°	@	Test	CREA	Method	Fixed Time	Direction	Ascend
Unit	mg/dL	Decimal	1	Prim. Wavelength	510	Sec. Wavelength	670
Sample Volume	20	R1 Volume	200	R2 Volume	0	Line Limit	20
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	2 -- 12	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0	Print Number	@
Upper	@	Upper	@	Upper	12.0		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	CREA	Method	Fixed Time	Direction	Ascend
Unit	mg/dL	Decimal	1	Prim. Wavelength	510	Sec. Wavelength	
Sample Volume	20	R1 Volume	160	R2 Volume	40	Line Limit	12
				Response			
Incubation	3	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	3 -- 8	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0	Print Number	@
Upper	@	Upper	@	Upper	12.0		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	CREA	Method	Fixed Time	Direction	Ascend
Unit	mg/dL	Decimal	1	Prim. Wavelength	510	Sec. Wavelength	670
Sample Volume	20	R1 Volume	200	R2 Volume	0	Line Limit	20
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	5 -- 15	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0	Print Number	@
Upper	@	Upper	@	Upper	10.0		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	FERR	Method	Fixed Time	Direction	Ascend
Unit	µg/L	Decimal	0	Prim. Wavelength	546	Sec. Wavelength	0
Sample Volume	6	R1 Volume	200	R2 Volume	0	Line Limit	20
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	3 -- 28	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0	Print Number	@
Upper	@	Upper	@	Upper	500		
Calibration							
Method	Curve	Std. 0	0				
		Std. 1	#				
		Std. 2	#				
		Std. 3	#				
		Std. 4	#				
		Std. 5	#				

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PARAMETERS							
N°	@	Test	Fe-Ferrozi	Method	End Point	Direction	Ascend
Unit	µg/dL	Decimal	1	Prim. Wavelength	578	Sec. Wavelength	0
Sample Volume	25	R1 Volume	200	R2 Volume	50	Line Limit	N/A
				Response			
Incubation	2	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	-2 - 25	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0	Print Number	@
Upper	@	Upper	@	Upper	1000		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	FALC	Method	Kinetic	Direction	Ascend
Unit	U/L	Decimal	0	Prim. Wavelength	405	Sec. Wavelength	0
Sample Volume	4	R1 Volume	200	R2 Volume	0	Line Limit	20
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	3 -- 10	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0	Print Number	@
Upper	@	Upper	@	Upper	1500		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	FOSF	Method	End Point	Direction	Ascend
Unit	mg/dL	Decimal	1	Prim. Wavelength	340	Sec. Wavelength	700
Sample Volume	3	R1 Volume	300	R2 Volume		Line Limit	N/A
				Response			
Incubation		Check Antigen	No	Lower	--	Sample Volume	@
Reaction	0 -- 25	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0.1	Print Number	@
Upper	@	Upper	@	Upper	20		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	FR	Method	End-Point	Direction	Ascend
Unit	UI/mL	Decimal	0	Prim. Wavelength	630	Sec. Wavelength	0
Sample Volume	3	R1 Volume	240	R2 Volume	60	Line Limit	N/A
				Response			
Incubation	5	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	-1 -- 10	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	2	Print Number	@
Upper	@	Upper	@	Upper	160		
Calibration							
Method	Curve	Std. 0	0				
		Std. 1	#				
		Std. 2	#				
		Std. 3	#				
		Std. 4	#				
		Std. 5	#				

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PARAMETERS							
N°	@	Test	FRUT	Method	Fixed Time	Direction	Ascend
Unit	mmol/L	Decimal	1	Prim. Wavelength	546	Sec. Wavelength	0
Sample Volume	10	R1 Volume	200	R2 Volume	0	Line Limit	20
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	50 -- 75	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0.1	Print Number	@
Upper	@	Upper	@	Upper	7.0		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	GLIC	Method	End Point	Direction	Ascend
Unit	mg/dL	Decimal	0	Prim. Wavelength	510	Sec. Wavelength	700
Sample Volume	3	R1 Volume	300	R2 Volume	0	Line Limit	N/A
Incubation	--	Check Antigen	No	Response		Sample Volume	@
Reaction	0 -- 50	Substrate	N/A	Lower	--	Dilution	@
R1 Blank		Mix R Blank		Upper	--		
Lower	@	Lower	@	Linearity		Full Name	@
Upper	@	Upper	@	Lower	1	Print Number	@
				Upper	500		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	HDL	Method	End Point	Direction	Ascend
Unit	mg/dL	Decimal	0	Prim. Wavelength	578	Sec. Wavelength	700
Sample Volume	3	R1 Volume	225	R2 Volume	75	Line Limit	N/A
				Response			
Incubation	25	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	-1 - 25	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0	Print Number	@
Upper	@	Upper	@	Upper	200		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	LDLD	Method	End Point	Direction	Ascend
Unit	mg/dL	Decimal	0	Prim. Wavelength	546	Sec. Wavelength	700
Sample Volume	2	R1 Volume	200	R2 Volume	67	Line Limit	N/A
				Response			
Incubation	25	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	-1 -- 25	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	1	Print Number	@
Upper	@	Upper	@	Upper	990		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	LDH	Method	Kinetic	Direction	Descend
Unit	U/L	Decimal	0	Prim. Wavelength	340	Sec. Wavelength	670
Sample Volume	4	R1 Volume	200	R2 Volume	0	Line Limit	20
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	3 -- 10	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0	Print Number	@
Upper	@	Upper	@	Upper	2000		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	Lipase	Method	Kinetic	Direction	Ascend
Unit	U/L	Decimal	0	Prim. Wavelength	578	Sec. Wavelength	700
Sample Volume	3	R1 Volume	200	R2 Volume	114	Line Limit	N/A
				Response			
Incubation	2	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	8 - 16	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0	Print Number	@
Upper	@	Upper	@	Upper	286		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	Mg	Method	End Point	Direction	Ascend
Unit	mg/dL	Decimal	1	Prim. Wavelength	510	Sec. Wavelength	0
Sample Volume	3	R1 Volume	300	R2 Volume	0	Line Limit	N/A
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	0 -- 10	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0.1	Print Number	@
Upper	@	Upper	@	Upper	4.5		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	Mau	Method	Fixed-Time	Direction	Ascend
Unit	mg/L	Decimal	0	Prim. Wavelength	546	Sec. Wavelength	0
Sample Volume	3	R1 Volume	420	R2 Volume	0	Line Limit	20
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	3 -- 14	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	1	Print Number	@
Upper	@	Upper	@	Upper	130		
Calibration							
Method	Linear	Std. 0	0.0 (NaCl)				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	Mau	Method	Fixed-Time	Direction	Ascend
Unit	mg/L	Decimal	0	Prim. Wavelength	546	Sec. Wavelength	0
Sample Volume	3	R1 Volume	340	R2 Volume	85	Line Limit	20
				Response			
Incubation	3	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	4 -- 14	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	1	Print Number	@
Upper	@	Upper	@	Upper	130		
Calibration							
Method	Linear	Std. 0	0.0 (NaCl)				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	PCR-AS	Method	Fixed-time	Direction	Ascend
Unit	mg/L	Decimal	2	Prim. Wavelength	546	Sec. Wavelength	0
Sample Volume	3	R1 Volume	225	R2 Volume	0	Line Limit	20
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	3 -- 28	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0.06	Print Number	@
Upper	@	Upper	@	Upper	15		
Calibration							
Method	Curve	Std. 0	0.0				
		Std. 1	#				
		Std. 2	#				
		Std. 3	#				
		Std. 4	#				
		Std. 5	#				

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PARAMETERS							
N°	@	Test	PROT	Method	End Point	Direction	Ascend
Unit	g/dL	Decimal	1	Prim. Wavelength	546	Sec. Wavelength	670
Sample Volume	4	R1 Volume	200	R2 Volume	0	Line Limit	N/A
Incubation	--	Check Antigen	No	Response	Lower	Sample Volume	@
Reaction	0 - 34	Substrate	N/A		Upper	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0	Print Number	@
Upper	@	Upper	@	Upper	14		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	Pt Urina	Method	End Point	Direction	Ascend
Unit	mg/dL	Decimal	0	Prim. Wavelength	578	Sec. Wavelength	670
Sample Volume	10	R1 Volume	200	R2 Volume	0	Line Limit	N/A
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	0 -- 25	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0	Print Number	@
Upper	@	Upper	@	Upper	100		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	TRI	Method	End Point	Direction	Ascend
Unit	mg/dL	Decimal	0	Prim. Wavelength	510	Sec. Wavelength	700
Sample Volume	3	R1 Volume	300	R2 Volume	0	Line Limit	N/A
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	0 -- 50	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	1	Print Number	@
Upper	@	Upper	@	Upper	1100		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	URE	Method	Fixed Time	Direction	Descend
Unit	mg/dL	Decimal	1	Prim. Wavelength	340	Sec. Wavelength	670
Sample Volume	3	R1 Volume	300	R2 Volume	0	Line Limit	20
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	2 - 6	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0	Print Number	@
Upper	@	Upper	@	Upper	300		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

REVISÃO: 03/11