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Instrument Settings

Test Methodology	
Kind of Process	End Point
Type	Linear
Filters	510 / 630
Reaction Direction	Increase
Reagent #1 (µL)	300
Reagent #2 (µL)	-
Sample Starter	Inactive
Delay Time (sec.)	0
Incubation Time (sec.)	280
Reading Time (sec.)	20

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Instrument Settings

Test Methodology	Turbidimetry
Kind of Process	Fixed time
Type	Linear
Filters	546
Reaction Direction	Increase
Reagent #1 (µL)	300
Reagent #2 (µL)	-
Sample Starter	Inactive
Delay Time (sec.)	10
Incubation Time (sec.)	30 - 0
Reading Time (sec.)	120

Units – Serum	UI/mL
Units – Urine	
Number of needle washes	2/2
Number of cuvette washes	2
Dynamic Blank	Inactive
Reagent Blank	Every Run
Reagent Limit (mABS)	1500
Curve Acceptance (%)	100
Instrument Factor	0.000
Shift	0.000

SERUM	
Name	AEO
Sample (µL)	3
Pre-Dilution	1.00
<u>Dilution</u>	
Factor	2.00
Test Limit	800
Initial ABS	0
Final ABS	1900
Max ABS Delta	500
<u>Normal Values</u>	
Man	0 – 200
Woman	0 – 200
Child	0 – 150
Calibrator	
Std. 1	0
Std. 2	#

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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Instrument Settings

Test Methodology	
Kind of Process	Kinetic
Type	Linear
Filters	340
Reaction Direction	Decrease
Reagent #1 (µL)	300
Reagent #2 (µL)	----
Sample Starter	Inactive
Delay Time (sec.)	0
Incubation Time (sec.)	60
Reading Time (sec.)	210

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Instrument Settings

Test Methodology	
Kind of Process	Kinetic
Type	Linear
Filters	405
Reaction Direction	Increase
Reagent #1 (µL)	250
Reagent #2 (µL)	----
Sample Starter	Inactive
Delay Time (sec.)	0
Incubation Time (sec.)	60
Reading Time (sec.)	240

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Instrument Settings

Test Methodology	DCA
Kind of Process	End Point 2 points
Type	Linear
Filters	546/630
Reaction Direction	Increase
Reagent #1 (µL)	300
Reagent #2 (µL)	75
Sample Starter	Inactive
Delay Time (sec.)	0
Incubation Time (sec.)	20 - 20
Reading Time (sec.)	180

Units – Serum	mg/dL
Units – Urine	---
Number of needle washes	1/1
Number of cuvette washes	1
Dynamic Blank	Inactive
Reagent Blank	Every Run
Reagent Limit (mABS)	300
Curve Acceptance (%)	20
Instrument Factor	0.000
Shift	0.000

SERUM	
Name	BIL-D
Sample (µL)	30
Pré-Dilution	1.00
<u>Dilution</u>	
Factor	10.00
Test Limit	12
Initial ABS	@
Final ABS	@
Max ABS Delta	99999
<u>Normal Values</u>	
Man	0 – 0.2
Woman	0 – 0.2
Child	-
Calibrator	
Std. 1	0
Std. 2	#

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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Instrument Settings

Test Methodology	
Kind of Process	End Point 2 points
Type	Linear
Filters	546/630
Reaction Direction	Increase
Reagent #1 (µL)	300
Reagent #2 (µL)	75
Sample Starter	Inactive
Delay Time (sec.)	0
Incubation Time (sec.)	30 - 240
Reading Time (sec.)	20

Units – Serum	mg/dL
Units – Urine	
Number of needle washes	1/1
Number of cuvette washes	1
Dynamic Blank	Inactive
Reagent Blank	Every Run
Reagent Limit (mABS)	300
Curve Acceptance (%)	20
Instrument Factor	0.000
Shift	0.000

SERUM	
Name	BIL-T
Sample (µL)	19
Pré-Dilution	1.00
<u>Dilution</u>	
Factor	10.00
Test Limit	30
Initial ABS	@
Final ABS	@
Max ABS Delta	99999
<u>Normal Values</u>	
Man	0 – 1.2
Woman	0 – 1.2
Child	-
Calibrator	
Std. 1	0
Std. 2	#

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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Instrument Settings

Test Methodology	
Kind of Process	End Point
Type	Linear
Filters	578
Reaction Direction	Increase
Reagent #1 (µL)	300
Reagent #2 (µL)	-
Sample Starter	Inactive
Delay Time (sec.)	0
Incubation Time (sec.)	30
Reading Time (sec.)	10

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Instrument Settings

Test Methodology	
Kind of Process	End Point
Type	Linear
Filters	510/630
Reaction Direction	Increase
Reagent #1 (µL)	300
Reagent #2 (µL)	-
Sample Starter	Inactive
Delay Time (sec.)	0
Incubation Time (sec.)	300
Reading Time (sec.)	20

Units – Serum	mg/dL
Units – Urine	
Number of needle washes	1/1
Number of cuvette washes	1
Dynamic Blank	Inactive
Reagent Blank	Every Run
Reagent Limit (mABS)	250
Curve Acceptance (%)	20
Instrument Factor	0.000
Shift	0.000

SERUM	
Name	CHOL
Sample (µL)	3
Pré-Dilution	1.00
<u>Dilution</u>	
Factor	2.00
Test Limit	500
Initial ABS	@
Final ABS	@
Max ABS Delta	99999
<u>Normal Values</u>	
Man	< 200
Woman	< 200
Child	-----
Calibrator	
Std. 1	0
Std. 2	#

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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Test Methodology	Tiocianato
Kind of Process	End Point
Type	Linear
Filters	480/700
Reaction Direction	Increase
Reagent #1 (µL)	300
Reagent #2 (µL)	-
Sample Starter	Inactive
Delay Time (sec.)	0
Incubation Time (sec.)	120
Reading Time (sec.)	10

Units – Serum	mEq/L
Units – Urine	mEq/L
Number of needle washes	2
Number of cuvette washes	2
Dynamic Blank	Inactive
Reagent Blank	Every Run
Reagent Limit (mABS)	1500
Curve Acceptance (%)	100
Instrument Factor	0.000
Shift	0.000

SERUM	
Name	Cloretos
Sample (µL)	3
Pré-Dilution	1.00
<u>Dilution</u>	
Factor	2.00
Test Limit	130
Initial ABS	
Final ABS	
Max ABS Delta	99999
Re-run Hyperactive	Active
Re-run Pathological	Inactive
<u>Normal Values</u>	
Man	98 – 110
Woman	98 – 110
Child	98 – 110
Calibrator	
Std. 1	(*)

URINE	
Name	Cloretos
Sample (µL)	3
Pré-Dilution	1.00
<u>Dilution</u>	
Factor	2.0
Test Limit	130
Initial ABS	
Final ABS	
Max ABS Delta	99999
Re-run Hyperactive	Active
Re-run Pathological	Inactive
<u>Normal Values</u>	
Man	
Woman	
Child	

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Instrument Settings

Test Methodology	
Kind of Process	Fixed Time
Type	Linear
Filters	510
Reaction Direction	Increase
Reagent #1 (µL)	60
Reagent #2 (µL)	240
Sample Starter	Inactive
Delay Time (sec.)	0
Incubation Time (sec.)	60 - 0
Reading Time (sec.)	120

Units – Serum	mg/dL
Units – Urine	
Number of needle washes	2/2
Number of cuvette washes	2
Dynamic Blank	Inactive
Reagent Blank	Every Run
Reagent Limit (mABS)	1500
Curve Acceptance (%)	100
Instrument Factor	0.000
Shift	0.000

SERUM	
Name	CREA
Sample (µL)	30
Pré-Dilution	1.00
<u>Dilution</u>	
Factor	2.00
Test Limit	12.0
Initial ABS	@
Final ABS	@
Max ABS Delta	500
<u>Normal Values</u>	
Man	0.6 – 1.3
Woman	0.6 – 1.1
Child	-----
Calibrator	
Std. 1	0
Std. 2	#

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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Instrument Settings

Test Methodology	
Kind of Process	Fixed Time
Type	Linear
Filters	510
Reaction Direction	Increase
Reagent #1 (µL)	300
Reagent #2 (µL)	-
Sample Starter	Inactive
Delay Time (sec.)	0
Incubation Time (sec.)	60
Reading Time (sec.)	120

SERUM	
Name	CREA
Sample (µL)	30
Pré-Dilution	1.00
<u>Dilution</u>	
Factor	2.00
Test Limit	10.0
Initial ABS	@
Final ABS	@
Max ABS Delta	500
<u>Normal Values</u>	
Man	0.9 – 1.3
Woman	0.6 – 1.1
Child	-----
Calibrator	
Std. 1	0
Std. 2	#

Units – Serum	mg/dL
Units – Urine	
Number of needle washes	2/2
Number of cuvette washes	2
Dynamic Blank	Inactive
Reagent Blank	Every Run
Reagent Limit (mABS)	1500
Curve Acceptance (%)	100
Instrument Factor	0.000
Shift	0.000

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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Reagente 1: Tampão (1)

Reagente 2: Látex Ferritina (2)

Test Methodology	Turbidimetry
Method	End Point
Kind of process	Cubic Spline
Filters	546
Reaction Direction	Increasing
Reagent #1 (µL)	200
Reagent #2 (µL)	100
Sample Starter	Inactive
Delay Time (sec.)	0
Incubation Time (sec.)	0/30
Reading Time (sec.)	270

Units – Serum	µg/L
Units – Urine	
Number of needle washes	2/2
Number of cuvette washes	2
Dynamic Blank	Inactive
Reagent Blank	Every Run
Reagent Limit (mABS)	1600
Curve Acceptance (%)	100
Instrument Factor	0.000
Shift	0.000

SERUM	
Name	FERRI
Sample (µL)	9
Pré-Dilution	1.00
<u>Dilution</u>	
Factor	2.00
Test Limit	500
Initial ABS	1900
Final ABS	1900
Max ABS Delta	500
<u>Normal Values</u>	
Man	20 – 250
Woman	20 – 200
Child	7 - 140
Calibrator	
Std. 1	0 µg/L
Std. 2	#
Std. 3	#
Std. 4	#
Std. 5	#
Std. 6	#

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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Instrument Settings

Test Methodology	
Kind of Process	End Point 2 points
Type	Linear
Filters	578
Reaction Direction	Increase
Reagent #1 (µL)	240
Reagent #2 (µL)	60
Sample Starter	Inactive
Delay Time (sec.)	0
Incubation Time (sec.)	60 - 0
Reading Time (sec.)	300

Units – Serum	µg/dL
Units – Urine	
Number of needle washes	2/2
Number of cuvette washes	2
Dynamic Blank	Inactive
Reagent Blank	Every Run
Reagent Limit (mABS)	1500
Curve Acceptance (%)	100
Instrument Factor	0.000
Shift	0.000

SERUM	
Name	Fe
Sample (µL)	30
Pré-Dilution	1.00
<u>Dilution</u>	
Factor	2.00
Test Limit	1000
Initial ABS	1900
Final ABS	1900
Max ABS Delta	500
<u>Normal Values</u>	
Man	65 - 170
Woman	50 - 170
Child	-----
Calibrator	
Std. 1	0
Std. 2	#

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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Instrument Settings

Test Methodology	
Kind of Process	Kinetc
Type	Linear
Filters	405
Reaction Direction	Increase
Reagent #1 (µL)	300
Reagent #2 (µL)	----
Sample Starter	Inactive
Delay Time (sec.)	0
Incubation Time (sec.)	60
Reading Time (sec.)	180

Units – Serum	U/L
Units – Urine	
Number of needle washes	2/2
Number of cuvette washes	2
Dynamic Blank	Inactive
Reagent Blank	Every Run
Reagent Limit (mABS)	1500
Curve Acceptance (%)	100
Instrument Factor	0.000
Shift	0.000

SERUM	
Name	FALC
Sample (µL)	6
Pré-Dilution	1.00
<u>Dilution</u>	
Factor	2.00
Test Limit	1500
Initial ABS	@
Final ABS	@
Max ABS Delta	500
<u>Normal Values</u>	
Man	27 – 100
Woman	27 – 100
Child	-
Calibrator	
Std. 1	0
Std. 2	#

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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Test Methodology	Molibd UV
Method	End Point
Kind of Process	Linear
Filters	340/700
Reaction Direction	Increase
Reagent #1 (µL)	300
Reagent #2 (µL)	0
Sample Starter	(Inactive)
Delay Time (sec.)	0
Incubation Time (sec.)	300
Reading Time (sec.)	10

Units – Serum	mg/dL
Units – Urine	mg/dL
Number of needle washes	2
Number of cuvette washes	2
Dynamic Blank	Inactive
Reagent Blank	Every Run
Reagent Limit (mABS)	1500
Curve Acceptance (%)	100

SERUM	
Name	Fósforo
Sample (µL)	3
Pré-Dilution	1.00
<u>Dilution</u>	
Factor	2.00
Test Limit (conc)	20
Max ABS Delta (mABS)	2000
Check Prozone	inactive
Instrumental Factor	0.000
Shift	0.000
Re-Run Hyperactive	active
Re-Run Pathological	Inactive
<u>Normal range</u>	
Male	2.50 – 4.80
Female	2.50 – 4.80
Child	4.50 – 6.70

URINE	
Name	Fósforo
Sample (µL)	3
Pré-Dilution	20.00
<u>Dilution</u>	
Factor	1.00
Test Limit (conc)	20
Check Prozone	inactive
Max ABS Delta (mABS)	2000
Instrumental Factor	0.000
Shift	0.000
Re-Run Hyperactive	active
Re-Run Pathological	Inactive
<u>Normal range</u>	
Male	0.00 / 0.00
Female	0.00 / 0.00
Child	0.00 / 0.00

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Instrument Settings

Test Methodology	
Kind of Process	Fixed Time
Type	Linear
Filters	546
Reaction Direction	Increase
Reagent #1 (µL)	300
Reagent #2 (µL)	----
Sample Starter	Inactive
Delay Time (sec.)	0
Incubation Time (sec.)	600
Reading Time (sec.)	300

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Instrument Settings

Test Methodology	
Kind of Process	End Point
Type	Linear
Filters	510
Reaction Direction	Increase
Reagent #1 (µL)	300
Reagent #2 (µL)	-
Sample Starter	Inactive
Delay Time (sec.)	0
Incubation Time (sec.)	300
Reading Time (sec.)	10

Units – Serum	mg/dL
Units – Urine	
Number of needle washes	1
Number of cuvette washes	1
Dynamic Blank	Inactive
Reagent Blank	Every Run
Reagent Limit (mABS)	250
Curve Acceptance (%)	20
Instrument Factor	0.000
Shift	0.000

SERUM	
Name	GLU
Sample (µL)	3
Pré-Dilution	1.00
<u>Dilution</u>	
Factor	2.00
Test Limit	500
Initial ABS	@
Final ABS	@
Max ABS Delta	500
<u>Normal Values</u>	
Man	70 - 105
Woman	70 - 105
Child	70 - 90
Calibrator	
Std. 1	0
Std. 2	#

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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Test Methodology	Turbidimetry
Kind of Process	End Point
Type	Linear
Filters	546/700
Reaction Direction	Increase
Reagent #1 (µL)	300
Reagent #2 (µL)	-
Sample Starter	Inactive
Delay Time (sec.)	10
Incubation Time (sec.)	30 - 0
Reading Time (sec.)	300

Units – Serum	mmol/L
Units – Urine	
Number of needle washes	2/2
Number of cuvette washes	2
Dynamic Blank	Inactive
Reagent Blank	Every Run
Reagent Limit (mABS)	1500
Curve Acceptance (%)	100
Instrument Factor	0.000
Shift	0.000

SERUM	
Name	Hb
Sample (µL)	50
Pré-Dilution	1.00
<i>Dilution</i>	
Factor	2.00
Test Limit	24
Initial ABS	0
Final ABS	1900
Max ABS Delta	500
<i>Normal Values</i>	
Man	
Woman	
Child	
Calibrator	
Std. 1	0.000
Std. 2	#

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<i>Dilution</i>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<i>Normal Values</i>	
Man	
Woman	
Child	

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Test Methodology	Turbidimetry
Kind of Process	End Point 2 points
Type	Cubic Spline
Filters	340
Reaction Direction	Increase
Reagent #1 (µL)	250
Reagent #2 (µL)	50
Sample Starter	Inactive
Delay Time (sec.)	10
Incubation Time (sec.)	300 - 0
Reading Time (sec.)	660

Units – Serum	mmol/L
Units – Urine	
Number of needle washes	2/2
Number of cuvette washes	2
Dynamic Blank	Inactive
Reagent Blank	Every Run
Reagent Limit (mABS)	1600
Curve Acceptance (%)	100
Instrument Factor	0.000
Shift	0.000

SERUM	
Name	HbA1c
Sample (µL)	10
Pré-Dilution	1.00
<i>Dilution</i>	
Factor	2.00
Test Limit	2.50
Initial ABS	1900
Final ABS	1900
Max ABS Delta	500
<i>Normal Values</i>	
Man	
Woman	
Child	
Calibrator	
Std. 1	0.000
Std. 2	#
Std. 3	#
Std. 4	#
Std.5	#
Std.6	#

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<i>Dilution</i>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<i>Normal Values</i>	
Man	
Woman	
Child	

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Instrument Settings

Test Methodology	
Kind of Process	End Point 2 points
Type	Linear
Filters	578/700
Reaction Direction	Increase
Reagent #1 (µL)	300
Reagent #2 (µL)	100
Sample Starter	Inactive
Delay Time (sec.)	0
Incubation Time (sec.)	300 - 0
Reading Time (sec.)	300

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Instrument Settings

Test Methodology	
Kind of Process	End Point 2 points
Type	Linear
Filters	546
Reaction Direction	Increase
Reagent #1 (µL)	300
Reagent #2 (µL)	100
Sample Starter	Inactive
Delay Time (sec.)	0
Incubation Time (sec.)	300 - 0
Reading Time (sec.)	300

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Test Methodology	Colorimetric
Method	Fixed Time
Kind of process	Linear
Filters	578
Reaction Direction	Increasing
Reagent #1 (µL)	200
Reagent #2 (µL)	114
Sample Starter	Inactive
Delay Time (sec.)	30
Incubation Time (sec.)	0 - 90
Reading Time (sec.)	90

Units – Serum	U/L
Units – Urine	
Number of needle washes	2/2
Number of cuvette washes	2
Dynamic Blank	Inactive
Reagent Blank	Every Run
Reagent Limit (mABS)	2000
Curve Acceptance (%)	100
Instrument Factor	0.000
Shift	0.000

SERUM	
Name	Lipase
Sample (µL)	3
Pré-Dilution	1;2
<u>Dilution</u>	
Factor	1.00
Test Limit	286
Initial ABS	@
Final ABS	@
Max ABS Delta	2000
<u>Normal Values</u>	
Man	13/60
Woman	13/60
Child	-----
Calibrator	
Std. 1	#

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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Instrument Settings

Test Methodology	
Kind of Process	End Point
Type	Linear
Filters	510
Reaction Direction	Increase
Reagent #1 (µL)	300
Reagent #2 (µL)	-
Sample Starter	Inactive
Delay Time (sec.)	10
Incubation Time (sec.)	120
Reading Time (sec.)	20

SERUM	
Name	MG
Sample (µL)	3
Pré-Dilution	1.00
<u>Dilution</u>	
Factor	2.00
Test Limit	4.5
Initial ABS	@
Final ABS	@
Max ABS Delta	99999
<u>Normal Values</u>	
Man	1.9 – 2.5
Woman	1.9 – 2.5
Child	@
Calibrator	
Std. 1	0
Std. 2	#

Units – Serum	mg/dL
Units – Urine	
Number of needle washes	2
Number of cuvette washes	2
Dynamic Blank	Inactive
Reagent Blank	Every Run
Reagent Limit (mABS)	800
Curve Acceptance (%)	100
Instrument Factor	0.000
Shift	0.000

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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Test Methodology	Turbidimetri
Kind of Process	Fixed Time
Type	Linear
Filters	546
Reaction Direction	Increase
Reagent #1 (µL)	320
Reagent #2 (µL)	80
Sample Starter	Inactive
Delay Time (sec.)	10
Incubation Time (sec.)	30 - 0
Reading Time (sec.)	120

SERUM	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	
Calibrator	
Std. 1	0
Std. 2	#

Units – Serum	
Units – Urine	mg/L
Number of needle washes	2/2
Number of cuvette washes	2
Dynamic Blank	Inactive
Reagent Blank	Every Run
Reagent Limit (mABS)	1500
Curve Acceptance (%)	100
Instrument Factor	0.000
Shift	0.000

URINE	
Name	MAU
Sample (µL)	3
Pré-Dilution	1.00
<u>Dilution</u>	
Factor	2.00
Test Limit	130
Initial ABS	1900
Final ABS	1900
Max ABS Delta	500
<u>Normal Values</u>	
Man	0 – 15
Woman	0 – 15
Child	

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Test Methodology	Turbidimetri
Kind of Process	Fixed Time
Type	Linear
Filters	546
Reaction Direction	Increase
Reagent #1 (µL)	420
Reagent #2 (µL)	
Sample Starter	Inactive
Delay Time (sec.)	0
Incubation Time (sec.)	0
Reading Time (sec.)	120

SERUM	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	
Calibrator	
Std. 1	0
Std. 2	#

Units – Serum	
Units – Urine	mg/L
Number of needle washes	2/2
Number of cuvette washes	2
Dynamic Blank	Inactive
Reagent Blank	Every Run
Reagent Limit (mABS)	1500
Curve Acceptance (%)	100
Instrument Factor	0.000
Shift	0.000

URINE	
Name	MAU
Sample (µL)	3
Pré-Dilution	1.00
<u>Dilution</u>	
Factor	2.00
Test Limit	130
Initial ABS	1900
Final ABS	1900
Max ABS Delta	500
<u>Normal Values</u>	
Man	0 – 15
Woman	0 – 15
Child	

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Instrument Settings

Test Methodology	Turbidimetry
Kind of Process	Fixed Time
Type	Linear
Filters	546/700
Reaction Direction	Increase
Reagent #1 (µL)	240
Reagent #2 (µL)	60
Sample Starter	Inactive
Delay Time (sec.)	30
Incubation Time (sec.)	0/10
Reading Time (sec.)	120

Units – Serum	mg/L
Units – Urine	
Number of needle washes	2/2
Number of cuvette washes	2
Dynamic Blank	Inactive
Reagent Blank	Every Run
Reagent Limit (mABS)	1500
Curve Acceptance (%)	100
Instrument Factor	0.000
Shift	0.000

SERUM	
Name	PCR
Sample (µL)	3
Pre-Dilution	1.00
<i>Dilution</i>	
Factor	2.00
Test Limit	100
Initial ABS	0
Final ABS	1900
Max ABS Delta	500
<i>Normal Values</i>	
Man	0 – 5
Woman	0 – 5
Child	0 – 5
Calibrator	
Std. 1	0
Std. 2	#

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<i>Dilution</i>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<i>Normal Values</i>	
Man	
Woman	
Child	

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Test Methodology	Turbidimetry
Method	Fixed Time
Kind of Process	Cubic spline
Filters	546
Reaction Direction	Increasing
Reagent #1 (µL)	240
Reagent #2 (µL)	60
Sample Starter	Inactive
Delay Time (sec.)	10
Incubation Time (sec.)	30/0
Reading Time (sec.)	300

Units – Serum	mg/L
Units – Urine	
Number of needle washes	2/2
Number of cuvette washes	2
Dynamic Blank	Inactive
Reagent Blank	Every day
Reagent Limit (mABS)	1500
Curve Acceptance (%)	100
Instrument Factor	0.000
Shift	0.000

SERUM	
Name	PCR AS
Sample (µL)	4
Pre-Dilution	1.00
<u>Dilution</u>	
Factor	2.00
Test Limit	
Initial ABS	0.000
Final ABS	1900
Max ABS Delta	500
<u>Normal Values</u>	
Man	
Woman	
Child	
Calibrator	
Std. 1	0.0 mg/L (sol. salina)
Std. 2	#
Std. 3	#
Std. 4	#
Std. 5	#
Std. 6	#

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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Instrument Settings

Test Methodology	
Kind of Process	End Point
Type	Linear
Filters	546
Reaction Direction	Increase
Reagent #1 (µL)	300
Reagent #2 (µL)	-
Sample Starter	Inactive
Delay Time (sec.)	0
Incubation Time (sec.)	600
Reading Time (sec.)	10

Units – Serum	g/dL
Units – Urine	
Number of needle washes	2/2
Number of cuvette washes	2
Dynamic Blank	Inactive
Reagent Blank	Every Run
Reagent Limit (mABS)	1500
Curve Acceptance (%)	100
Instrument Factor	0.000
Shift	0.000

SERUM	
Name	PROT
Sample (µL)	6
Pré-Dilution	1.00
<u>Dilution</u>	
Factor	2.00
Test Limit	14.0
Initial ABS	@
Final ABS	@
Max ABS Delta	500
<u>Normal Values</u>	
Man	6.0 – 8.0
Woman	6.0 – 8.0
Child	6.0 – 8.0
Calibrator	
Std. 1	0
Std. 2	#

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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Instrument Settings

Test Methodology	
Kind of Process	End Point
Type	Linear
Filters	578
Reaction Direction	Increase
Reagent #1 (µL)	300
Reagent #2 (µL)	-
Sample Starter	Inactive
Delay Time (sec.)	0
Incubation Time (sec.)	600
Reading Time (sec.)	10

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Instrument Settings

Test Methodology	
Kind of Process	End Point
Type	Linear
Filters	510 / 630
Reaction Direction	Increase
Reagent #1 (µL)	300
Reagent #2 (µL)	-
Sample Starter	Inactive
Delay Time (sec.)	0
Incubation Time (sec.)	300
Reading Time (sec.)	20

Units – Serum	mg/dL
Units – Urine	
Number of needle washes	1
Number of cuvette washes	1
Dynamic Blank	Inactive
Reagent Blank	Every Run
Reagent Limit (mABS)	500
Curve Acceptance (%)	20
Instrument Factor	0.000
Shift	0.000

SERUM	
Name	TRIG
Sample (µL)	3
Pré-Dilution	1.00
<u>Dilution</u>	
Factor	5.00
Test Limit	1100
Initial ABS	@
Final ABS	@
Max ABS Delta	99999
<u>Normal Values</u>	
Man	60 - 150
Woman	60 - 150
Child	60 - 150
Calibrator	
Std. 1	0
Std. 2	#

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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Instrument Settings

Test Methodology	
Kind of Process	End Point
Type	Linear
Filters	578
Reaction Direction	Increase
Reagent #1 (µL)	300
Reagent #2 (µL)	300
Sample Starter	Inactive
Delay Time (sec.)	0
Incubation Time (sec.)	300 / 10
Reading Time (sec.)	300

Units – Serum	mg/dL
Units – Urine	
Number of needle washes	3/3
Number of cuvette washes	3
Dynamic Blank	Inactive
Reagent Blank	Every Run
Reagent Limit (mABS)	500
Curve Acceptance (%)	20
Instrument Factor	0.000
Shift	0.000

SERUM	
Name	UREA UV
Sample (µL)	3
Pré-Dilution	1.00
<i>Dilution</i>	
Factor	1.00
Test Limit	300
Initial ABS	@
Final ABS	@
Max ABS Delta	99999
<i>Normal Values</i>	
Man	15 – 40
Woman	15 – 40
Child	-
Calibrator	
Std. 1	0
Std. 2	#

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<i>Dilution</i>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<i>Normal Values</i>	
Man	
Woman	
Child	

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Instrument Settings

Test Methodology	
Kind of Process	Fixed time
Type	Linear
Filters	340 / 630
Reaction Direction	Decrease
Reagent #1 (µL)	240
Reagent #2 (µL)	60
Sample Starter	Active
Delay Time (sec.)	0 / 120
Incubation Time (sec.)	60
Reading Time (sec.)	60

Units – Serum	mg/dL
Units – Urine	
Number of needle washes	1/1
Number of cuvette washes	1
Dynamic Blank	Inactive
Reagent Blank	Every Run
Reagent Limit (mABS)	900
Curve Acceptance (%)	20
Instrument Factor	0.000
Shift	0.000

SERUM	
Name	UREA UV
Sample (µL)	3
Pré-Dilution	1.00
<u>Dilution</u>	
Factor	10.00
Test Limit	300
Initial ABS	@
Final ABS	@
Max ABS Delta	1000
<u>Normal Values</u>	
Man	15 – 40
Woman	15 – 40
Child	-
Calibrator	
Std. 1	0
Std. 2	#

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

EDIÇÃO: 09/11

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© BT 3000 PLUS é marca registrada de seus proprietários.