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GENERAL	Test name	Ácido Úrico
	Analysis mode	Endpoint mon.
	Sample type	Serum
	Units	mg/dL
	Reaction type	Increasing
	Decimals	2
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Bichrom.
	Volumes	
	Sample	8 µL
	Reagent 1	300 µL
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	505
	Reference	670
Times	Reading 1	312 s
	Reading 2	-
	Reagent 2	-
	Postdilution factor	2
CALIBRATION	Type of calibration	Múltiplo (*)
	Calibrator replicates	2
	Blank replicates	1
	Calibration curve	-
OPTIONS	Blank absorbance limit	0.180
	Kinetic blank limit	-
	Linerity limit	25.00 mg/dL

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GENERAL	Test name	AEO Turbidim
	Analysis mode	Endpoint mon.
	Sample type	Serum
	Units	UI/mL
	Reaction type	Increasing
	Turbidim. Test	yes
	Decimals	0
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	bichrom.
	Volumes	
	Sample	3
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	535
	Reference	670
	Times	
	Reading 1	192 s
Reading 2	-	
Reagent 2	-	
Postdilution factor	2	
CALIBRATION	Type of calibration	specific (*)
	Calibrator number	1
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	0.700
	Kinetic blank limit	-
	Linerity limit	800

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GENERAL	Test name	Albumina
	Analysis mode	Endpoint mon.
	Sample type	Serum
	Units	g/dL
	Reaction type	Increasing
	Decimals	2
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	monochrom.
	Volumes	
	Sample	3
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	635
	Reference	
	Times	
	Reading 1	120 s
Reading 2	-	
Reagent 2	-	
Postdilution factor	2	
CALIBRATION	Type of calibration	Múltiplo (*)
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	0.200
	Kinetic blank limit	-
	Linerity limit	6.0

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GENERAL	Test name	ALT
	Analysis mode	Kinetic mon.
	Sample type	Serum
	Units	U/L
	Reaction type	decreasing
	Decimals	0
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monoch.
	Volumes	
	Sample	15
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	340
	Reference	-
	Times	
	Reading 1	96 s
Reading 2	264	
Reagent 2	-	
Postdilution factor	2	
CALIBRATION	Type of calibration	Múltiplo (*)
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	1.100
	Kinetic blank limit	-
	Linerity limit	400

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GENERAL	Test name	Amilase
	Analysis mode	Kinetic mon.
	Sample type	Serum
	Units	U/L
	Reaction type	Increasing
	Decimals	0
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monoch.
	Volumes	
	Sample	6
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	405
	Reference	-
	Times	
	Reading 1	72 s
Reading 2	216 s	
Reagent 2	-	
Postdilution factor	2	
CALIBRATION	Type of calibration	Múltiplo (*)
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	0.700
	Kinetic blank limit	-
	Linerity limit	2000

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GENERAL	Test name	AST
	Analysis mode	kinetic mon.
	Sample type	Serum
	Units	U/L
	Reaction type	decreasing
	Decimals	0
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monoch.
	Volumes	
	Sample	15
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	340
	Reference	-
	Times	
	Reading 1	96 s
Reading 2	264 s	
Reagent 2	-	
Postdilution factor	2	
CALIBRATION	Type of calibration	Múltiplo (*)
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	1.100
	Kinetic blank limit	-
	Linerity limit	400

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GENERAL	Test name	Cálcio
	Analysis mode	Endpoint mon.
	Sample type	Serum
	Units	mg/dL
	Reaction type	Increasing
	Decimals	1
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monoch.
	Volumes	
	Sample	3
	Reagent 1	300
	Reagent 2	----
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	560
	Reference	-
	Times	
	Reading 1	120 s
Reading 2	----	
Reagent 2	----	
Postdilution factor	2	
CALIBRATION	Type of calibration	Múltiplo (*)
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	2.000
	Kinetic blank limit	-
	Linerity limit	16.0

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GENERAL	Test name	Cálcio Arsenazo	
	Analysis mode	Endpoint mon.	
	Sample type	Serum	
	Units	mg/dL	
	Reaction type	Increasing	
	Decimals	1	
	Replicates	1	
	Name of assoc. constituent	-	
PROCEDURE	Volumes	Type of reading	Monoch.
		Sample	3
		Reagent 1	300
		Reagent 2	-
		Washing	1.2
	Filters	Predilution factor	-
		Main	670
		Reference	-
	Times	Reading 1	120 s
		Reading 2	-
		Reagent 2	-
		Postdilution factor	2
	CALIBRATION	Type of calibration	Múltiplo (*)
Calibrator replicates		3	
Blank replicates		3	
Calibration curve		-	
OPTIONS	Blank absorbance limit	1.000	
	Kinetic blank limit	-	
	Linerity limit	17	

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GENERAL	Test name	CLF
	Analysis mode	Differential Bir.
	Sample type	Serum
	Units	µg/dL
	Reaction type	decreasing
	Decimals	1
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monoch.
	Volumes	
	Sample	30
	Reagent 1	235
	Reagent 2	65
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	560
	Reference	-
	Times	
	Reading 1	288 s
	Reading 2	600 s
Reagent 2	312 s	
Postdilution factor	2	
CALIBRATION	Type of calibration	specific (*)
	Calibrator replicates	2
	Blank replicates	2
	Calibration curve	-
OPTIONS	Blank absorbance limit	0.350
	Kinetic blank limit	-
	Linerity limit	500

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GENERAL	Test name	Cloretos	
	Analysis mode	Endpoint mon.	
	Sample type	Serum	
	Units	mg/dL	
	Reaction type	Increasing	
	Decimals	0	
	Replicates	1	
	Name of assoc. constituent	-	
PROCEDURE	Volumes	Type of reading	Bichrom.
		Sample	3 µL
		Reagent 1	300 µL
		Reagent 2	-
		Washing	1.2
		Predilution factor	-
	Filters	Main	505
		Reference	670
	Times	Reading 1	120 s
		Reading 2	-
		Reagent 2	-
		Postdilution factor	-
CALIBRATION	Type of calibration	Múltiplo (*)	
	Calibrator replicates	1	
	Blank replicates	1	
	Calibration curve	-	
OPTIONS	Blank absorbance limit	0.200	
	Kinetic blank limit	-	
	Linerity limit	130 mg/dL	

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GENERAL	Test name	CKMB
	Analysis mode	Kinetic mon.
	Sample type	Serum
	Units	U/L
	Reaction type	Increasing
	Decimals	0
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monoch.
	Volumes	
	Sample	15
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	340
	Reference	-
	Times	
	Reading 1	312 s
Reading 2	600 s	
	Reagent 2	-
	Postdilution factor	-
CALIBRATION	Type of calibration	Specific (*)
	Factor	
	Calibrator replicates	
	Blank replicates	3
OPTIONS	Calibration curve	-
	Blank absorbance limit	0.500
	Kinetic blank limit	-
	Linerity limit	500

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GENERAL	Test name	CK NAC
	Analysis mode	Kinetic mon.
	Sample type	Serum
	Units	U/L
	Reaction type	Increasing
	Decimals	0
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monoch.
	Volumes	
	Sample	6
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	340
	Times	
	Reference	-
	Reading 1	120 s
Reading 2	240 s	
Reagent 2	-	
Postdilution factor	2	
CALIBRATION	Type of calibration	specific (*)
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	0.500
	Kinetic blank limit	-
	Linerity limit	2000

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GENERAL	Test name	CK NAC
	Analysis mode	Kinetic mon.
	Sample type	Serum
	Units	U/L
	Reaction type	Increasing
	Decimals	0
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monoch.
	Volumes	
	Sample	6
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	340
	Reference	-
	Times	
	Reading 1	120 s
Reading 2	240 s	
	Reagent 2	-
	Postdilution factor	2
CALIBRATION	Type of calibration	specific (*)
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	0.500
	Kinetic blank limit	-
	Linerity limit	2000

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GENERAL	Test name	Colesterol	
	Analysis mode	Endpoint mon.	
	Sample type	Serum	
	Units	mg/dL	
	Reaction type	Increasing	
	Decimals	0	
	Replicates	1	
	Name of assoc. constituent	-	
PROCEDURE	Volumes	Type of reading	Bichrom.
		Sample	3 µL
		Reagent 1	300 µL
		Reagent 2	-
		Washing	1.2
	Filters	Predilution factor	-
		Main	505
	Times	Reference	670
		Reading 1	600 s
		Reading 2	-
		Reagent 2	-
		Postdilution factor	-
CALIBRATION	Type of calibration	Múltiplo (*)	
	Calibrator replicates	1	
	Blank replicates	1	
	Calibration curve	-	
OPTIONS	Blank absorbance limit	0.180	
	Kinetic blank limit	-	
	Linerity limit	500 mg/dL	

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GENERAL	Test name	Colesterol HDL
	Analysis mode	Endpoint mon.
	Sample type	Serum
	Units	mg/dL
	Reaction type	Increasing
	Decimals	0
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Bichrom.
	Volumes	
	Sample	30 µL
	Reagent 1	300 µL
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	505
	Reference	670
	Times	
	Reading 1	600 s
Reading 2	-	
Reagent 2	-	
Postdilution factor	2	
CALIBRATION	Type of calibration	specific (*)
	Calibrator replicates	1
	Blank replicates	1
	Calibration curve	-
OPTIONS	Blank absorbance limit	0.180
	Kinetic blank limit	-
	Linerity limit	200 mg/dL

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GENERAL	Test name	CHE
	Analysis mode	Kinetic bir.
	Sample type	Serum
	Units	U/L
	Reaction type	Decreasing
	Decimals	0
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monoch.
	Volumes	
	Sample	4
	Reagent 1	200
	Reagent 2	50
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	405
	Reference	-
	Times	
	Reading 1	312 s
Reading 2	408 s	
	Reagent 2	192 s
	Postdilution factor	2
CALIBRATION	Type of calibration	Multiple (*)
	Factor	
	Calibrator replicates	3
	Blank replicates	3
OPTIONS	Calibration curve	-
	Blank absorbance limit	1.100
	Kinetic blank limit	-
	Linerity limit	20000

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GENERAL	Test name	Creatinina	
	Analysis mode	Fixed-time mon.	
	Sample type	Serum	
	Units	mg/dL	
	Reaction type	Increasing	
	Decimals	2	
	Replicates	1	
	Name of assoc. constituent	-	
PROCEDURE	Type of reading	Monoch.	
	Volumes	Sample	30 µL
		Reagent 1	300 µL
		Reagent 2	-
		Washing	1.2
		Predilution factor	-
	Filters	Main	505
		Reference	-
	Times	Reading 1	48 s
		Reading 2	120 s
		Reagent 2	-
	Postdilution factor	2	
CALIBRATION	Type of calibration	Múltiplo (*)	
	Calibrator replicates	3	
	Blank replicates	3	
	Calibration curve	-	
OPTIONS	Blank absorbance limit	0.200	
	Kinetic blank limit	-	
	Linerity limit	12.0 mg/dL	

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GENERAL	Test name	Creatinina
	Analysis mode	Fixed-time mon.
	Sample type	Serum
	Units	mg/dL
	Reaction type	Increasing
	Decimals	2
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monoch.
	Volumes	
	Sample	30 µL
	Reagent 1	300 µL
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	505
	Reference	-
	Times	
	Reading 1	48 s
Reading 2	120 s	
	Reagent 2	-
	Postdilution factor	2
CALIBRATION	Type of calibration	Múltiplo (*)
	Calibrator replicates	2
	Blank replicates	1
	Calibration curve	-
OPTIONS	Blank absorbance limit	0.200
	Kinetic blank limit	-
	Linerity limit	12 mg/dL

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GENERAL	Test name	Creatinina
	Analysis mode	Fixed-time mon.
	Sample type	Serum
	Units	mg/dL
	Reaction type	Increasing
	Decimals	2
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monoch.
	Volumes	
	Sample	30 µL
	Reagent 1	300 µL
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	505
	Reference	-
	Times	
	Reading 1	72 s
Reading 2	192 s	
	Reagent 2	-
	Postdilution factor	2
CALIBRATION	Type of calibration	Múltiplo (*)
	Calibrator replicates	2
	Blank replicates	1
	Calibration curve	-
OPTIONS	Blank absorbance limit	0.400
	Kinetic blank limit	-
	Linerity limit	10 mg/dL

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GENERAL	Test name	Ferritina
	Analysis mode	Fixed Time mon.
	Sample type	Serum
	Units	µg/L
	Reaction type	Increasing
	Turbidim. Test	yes
	Decimals	0
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	monoch.
	Volumes	
	Sample	8
	Reagent 1	240
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	535
	Reference	-
	Times	
	Reading 1	48 s
Reading 2	360 s	
Reagent 2	-	
Postdilution factor	-	
CALIBRATION	Type of calibration	specific (*)
	Calibrator number	5
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	Increasing Polygonal
OPTIONS	Blank absorbance limit	1.600
	Kinetic blank limit	-
	Linerity limit	-

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GENERAL	Test name	Ferro Ferrozina
	Analysis mode	Differential bir.
	Sample type	Serum
	Units	µg/dL
	Reaction type	Increasing
	Decimals	0
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monoch.
	Volumes	
	Sample	40 µL
	Reagent 1	240 µL
	Reagent 2	60 µL
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	560
	Reference	-
	Times	
	Reading 1	72 s
Reading 2	408 s	
Reagent 2	96 s	
Postdilution factor	2	
CALIBRATION	Type of calibration	Específico (*)
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	0.300
	Kinetic blank limit	-
	Linerity limit	750 µg/dL

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GENERAL	Test name	Fosfatase Alcalina AMP
	Analysis mode	Kinetic mon.
	Sample type	Serum
	Units	U/L
	Reaction type	Increasing
	Decimals	0
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monoch.
	Volumes	
	Sample	6
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	405
	Reference	-
	Times	
Reading 1	72 s	
Reading 2	192 s	
	Reagent 2	-
	Postdilution factor	2
CALIBRATION	Type of calibration	Múltiplo (*)
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	1.000
	Kinetic blank limit	-
	Linerity limit	1500

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GENERAL	Test name	Fósforo UV	
	Analysis mode	Endpoint mon.	
	Sample type	Serum	
	Units	mg/dL	
	Reaction type	Increasing	
	Decimals	2	
	Replicates	1	
	Name of assoc. constituent	-	
PROCEDURE	Volumes	Type of reading	Monoch.
		Sample	3
		Reagent 1	300
		Reagent 2	-
		Washing	1.2
	Filters	Predilution factor	-
		Main	340
		Reference	-
	Times	Reading 1	312 s
		Reading 2	-
		Reagent 2	-
		Postdilution factor	2
CALIBRATION	Type of calibration	Múltiplo (*)	
	Calibrator replicates	1	
	Blank replicates	1	
	Calibration curve	-	
OPTIONS	Blank absorbance limit	0.500	
	Kinetic blank limit	-	
	Linerity limit	20.0 mg/dL	

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GENERAL	Test name	FR Turbidim
	Analysis mode	Endpoint bir.
	Sample type	Serum
	Units	UI/mL
	Reaction type	Increasing
	Turbidim. Test	yes
	Decimals	0
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	monoch.
	Volumes	
	Sample	3
	Reagent 1	240
	Reagent 2	60
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	635
	Reference	-
Times		
Reading 1	144 s	
Reading 2	-	
Reagent 2	24 s	
Postdilution factor	2	
CALIBRATION	Type of calibration	specific (*)
	Calibrator number	5
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	Increasing Polygonal
OPTIONS	Blank absorbance limit	1.400
	Kinetic blank limit	-
	Linerity limit	-

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GENERAL	Test name	Frutosamina
	Analysis mode	Fixed-time mon.
	Sample type	Serum
	Units	mmol/L
	Reaction type	Increasing
	Decimals	2
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monoch.
	Volumes	
	Sample	15
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	535
	Times	
	Reference	-
	Reading 1	408 s
Reading 2	600 s	
Reagent 2	-	
Postdilution factor	2	
CALIBRATION	Type of calibration	Specific (*)
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	0.065
	Kinetic blank limit	-
	Linerity limit	7.0

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GENERAL	Test name	Gama GT
	Analysis mode	Kinetic mon.
	Sample type	Serum
	Units	U/L
	Reaction type	Increasing
	Decimals	0
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monoch.
	Volumes	
	Sample	15
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	405
	Reference	-
	Times	
	Reading 1	72 s
Reading 2	192 s	
Reagent 2	-	
Postdilution factor	2	
CALIBRATION	Type of calibration	Múltiplo (*)
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	1.450
	Kinetic blank limit	-
	Linerity limit	700

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GENERAL	Test name	Glicose
	Analysis mode	Endpoint mon.
	Sample type	Serum
	Units	mg/dL
	Reaction type	Increasing
	Decimals	0
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	bichrom.
	Volumes	
	Sample	3
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	505
	Reference	670
	Times	
	Reading 1	312 s
Reading 2	-	
Reagent 2	-	
Postdilution factor	2	
CALIBRATION	Type of calibration	Múltiplo (*)
	Calibrator replicates	1
	Blank replicates	1
	Calibration curve	-
OPTIONS	Blank absorbance limit	0.180
	Kinetic blank limit	-
	Linerity limit	500 mg/dL

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GENERAL	Test name	HDL Direto
	Analysis mode	Diferencial birreagente
	Sample type	Serum
	Units	mg/dL
	Reaction type	Increasing
	Decimals	0
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monocromática
	Volumes	
	Sample	4
	Reagent 1	300
	Reagent 2	100
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	560
	Reference	
	Times	
Reading 1	168 s	
Reading 2	480 s	
Reagent 2	192 s	
Postdilution factor	2	
CALIBRATION	Type of calibration	Specific (*)
	Calibrator replicates	1
	Blank replicates	2
	Calibration curve	-
OPTIONS	Blank absorbance limit	0.200
	Kinetic blank limit	-
	Linerity limit	150 mg/dL

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GENERAL	Test name	LDL Direto
	Analysis mode	Differential bir.
	Sample type	Serum
	Units	mg/dL
	Reaction type	Increasing
	Decimals	0
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monoch.
	Volumes	
	Sample	3
	Reagent 1	300
	Reagent 2	100
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	535
	Reference	-
	Times	
	Reading 1	288 s
Reading 2	480 s	
Reagent 2	312 s	
Postdilution factor	2	
CALIBRATION	Type of calibration	specific (*)
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	0.100
	Kinetic blank limit	-
	Linerity limit	990

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GENERAL	Test name	LDH
	Analysis mode	Kinetic mon.
	Sample type	Serum
	Units	U/L
	Reaction type	decreasing
	Decimals	0
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monoch.
	Volumes	
	Sample	6
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	340
	Reference	-
	Times	
	Reading 1	72 s
Reading 2	192 s	
	Reagent 2	-
	Postdilution factor	2
CALIBRATION	Type of calibration	Múltiplo (*)
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	0.800
	Kinetic blank limit	-
	Linerity limit	1500

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GENERAL	Test name	LIPASE
	Analysis mode	Fixed time Bir.
	Sample type	Serum
	Units	µg/dL
	Reaction type	Increasing
	Decimals	1
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monoch.
	Volumes	
	Sample	3
	Reagent 1	200
	Reagent 2	110
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	560
	Reference	-
	Times	
	Reading 1	120 s
Reading 2	216 s	
Reagent 2	24 s	
Postdilution factor	2	
CALIBRATION	Type of calibration	multiple (*)
	Calibrator replicates	2
	Blank replicates	2
	Calibration curve	-
OPTIONS	Blank absorbance limit	0.200
	Kinetic blank limit	-
	Linerity limit	300

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GENERAL	Test name	Magnésio
	Analysis mode	Endpoint mon.
	Sample type	Serum
	Units	mg/dL
	Reaction type	Increasing
	Decimals	2
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monoch.
	Volumes	
	Sample	3
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	505
	Reference	-
	Times	
	Reading 1	120 s
Reading 2	-	
Reagent 2	-	
Postdilution factor	2	
CALIBRATION	Type of calibration	Múltiplo (*)
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	0.600
	Kinetic blank limit	-
	Linerity limit	4.5

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GENERAL	Test name	Magnésio
	Analysis mode	Endpoint mon.
	Sample type	Serum
	Units	mg/dL
	Reaction type	Increasing
	Decimals	2
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monoch.
	Volumes	
	Sample	3
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	505
	Reference	-
	Times	
	Reading 1	600 s
Reading 2	-	
Reagent 2	-	
Postdilution factor	2	
CALIBRATION	Type of calibration	Múltiplo (*)
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	1.000
	Kinetic blank limit	-
	Linerity limit	3.5

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GENERAL	Test name	Microalb
	Analysis mode	Fixed-Time
	Sample type	Urine
	Units	mg/L
	Reaction type	Increasing
	Turbidim. Test	yes
	Decimals	0
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monoch.
	Volumes	
	Sample	3
	Reagent 1	420
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	535
	Times	
	Reference	-
	Reading 1	48 s
Reading 2	168 s	
Reagent 2	-	
Postdilution factor	1.4	
CALIBRATION	Type of calibration	specific (*)
	Calibrator number	1
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	1.200
	Kinetic blank limit	-
	Linerity limit	130

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GENERAL	Test name	PCR Turbid.
	Analysis mode	Endpoint mon.
	Sample type	Serum
	Units	mg/L
	Reaction type	Increasing
	Turbidim. Type	yes
	Decimals	0
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	bichrom.
	Volumes	
	Sample	3
	Reagent 1	440
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	535
	Reference	670
	Times	
	Reading 1	192 s
Reading 2	-	
Reagent 2	-	
Postdilution factor	1.4	
CALIBRATION	Type of calibration	specific (*)
	Calibrator number	1
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	0.700
	Kinetic blank limit	-
	Linerity limit	150

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GENERAL	Test name	Proteínas totais
	Analysis mode	Endpoint mon.
	Sample type	Serum
	Units	g/dL
	Reaction type	Increasing
	Decimals	1
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	bichrom.
	Volumes	
	Sample	6
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	535
	Reference	670
	Times	
	Reading 1	600 s
Reading 2	-	
Reagent 2	-	
Postdilution factor	2	
CALIBRATION	Type of calibration	Múltiplo (*)
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	0.200
	Kinetic blank limit	-
	Linerity limit	14

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GENERAL	Test name	Proteinúria
	Analysis mode	Endpoint mon.
	Sample type	urine
	Units	mg/dL
	Reaction type	Increasing
	Decimals	0
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monochrom.
	Volumes	
	Sample	15
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	600
	Reference	
	Times	
	Reading 1	312 s
Reading 2	-	
Reagent 2	-	
Postdilution factor	2	
CALIBRATION	Type of calibration	specific (*)
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	0.300
	Kinetic blank limit	-
	Linerity limit	100

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GENERAL	Test name	Triglicérides	
	Analysis mode	Endpoint mon.	
	Sample type	Serum	
	Units	mg/dL	
	Reaction type	Increasing	
	Decimals	0	
	Replicates	1	
	Name of assoc. constituent	-	
PROCEDURE	Type of reading	bichrom.	
	Volumes	Sample	3
		Reagent 1	300
		Reagent 2	-
		Washing	1.2
		Predilution factor	-
	Filters	Main	505
		Reference	670
	Times	Reading 1	600 s
		Reading 2	-
		Reagent 2	-
		Postdilution factor	2
CALIBRATION	Type of calibration	Múltiplo (*)	
	Calibrator replicates	1	
	Blank replicates	1	
	Calibration curve	-	
OPTIONS	Blank absorbance limit	0.180	
	Kinetic blank limit	-	
	Linerity limit	1100 mg/dL	

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GENERAL	Test name	Uréia UV
	Analysis mode	Fixed-time mon.
	Sample type	Serum
	Units	mg/dL
	Reaction type	decreasing
	Decimals	0
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monoch.
	Volumes	
	Sample	3
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	340
	Reference	-
	Times	
	Reading 1	48 s
Reading 2	120 s	
	Reagent 2	-
	Postdilution factor	2
CALIBRATION	Type of calibration	Múltiplo (*)
	Calibrator replicates	2
	Blank replicates	1
	Calibration curve	-
OPTIONS	Blank absorbance limit	0.600
	Kinetic blank limit	-
	Linerity limit	300

REVISÃO: 02/11

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© A15 é marca registrada de seus proprietários.