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# Analisa

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<b>GENERAL</b>	<b>Test name</b>	<b>Ácido Úrico</b>
	<b>Analysis mode</b>	<b>Endpoint mon.</b>
	<b>Sample type</b>	<b>Serum</b>
	<b>Units</b>	<b>mg/dL</b>
	<b>Reaction type</b>	<b>Increasing</b>
	<b>Decimals</b>	<b>2</b>
	<b>Replicates</b>	<b>1</b>
	<b>Name of assoc. constituent</b>	<b>-</b>
<b>PROCEDURE</b>	<b>Type of reading</b>	<b>Bichrom.</b>
	<b>Volumes</b>	
	<b>Sample</b>	<b>6 µL</b>
	<b>Reagent 1</b>	<b>300 µL</b>
	<b>Reagent 2</b>	<b>-</b>
	<b>Washing</b>	<b>1.2</b>
	<b>Predilution factor</b>	<b>-</b>
	<b>Filters</b>	
	<b>Main</b>	<b>505</b>
	<b>Reference</b>	<b>670</b>
	<b>Times</b>	
	<b>Reading 1</b>	<b>600 s</b>
<b>Reading 2</b>	<b>-</b>	
<b>Reagent 2</b>	<b>-</b>	
<b>Postdilution factor</b>	<b>2</b>	
<b>CALIBRATION</b>	<b>Type of calibration</b>	<b>Múltiplo (*)</b>
	<b>Calibrator replicates</b>	<b>2</b>
	<b>Blank replicates</b>	<b>1</b>
	<b>Calibration curve</b>	<b>-</b>
<b>OPTIONS</b>	<b>Blank absorbance limit</b>	<b>0.180</b>
	<b>Kinetic blank limit</b>	<b>-</b>
	<b>Linerity limit</b>	<b>20.00 mg/dL</b>

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<b>GENERAL</b>	<b>Test name</b>	<b>AEO Turbidim</b>
	<b>Analysis mode</b>	<b>Endpoint mon.</b>
	<b>Sample type</b>	<b>Serum</b>
	<b>Units</b>	<b>UI/mL</b>
	<b>Reaction type</b>	<b>Increasing</b>
	<b>Turbidim. Test</b>	<b>yes</b>
	<b>Decimals</b>	<b>0</b>
	<b>Replicates</b>	<b>1</b>
	<b>Name of assoc. constituent</b>	<b>-</b>
<b>PROCEDURE</b>	<b>Type of reading</b>	<b>bichrom.</b>
	<b>Volumes</b>	
	<b>Sample</b>	<b>3</b>
	<b>Reagent 1</b>	<b>300</b>
	<b>Reagent 2</b>	<b>-</b>
	<b>Washing</b>	<b>1.2</b>
	<b>Predilution factor</b>	<b>-</b>
	<b>Filters</b>	
	<b>Main</b>	<b>535</b>
	<b>Reference</b>	<b>670</b>
	<b>Times</b>	
	<b>Reading 1</b>	<b>192 s</b>
<b>Reading 2</b>	<b>-</b>	
<b>Reagent 2</b>	<b>-</b>	
<b>Postdilution factor</b>	<b>2</b>	
<b>CALIBRATION</b>	<b>Type of calibration</b>	<b>specific (*)</b>
	<b>Calibrator number</b>	<b>1</b>
	<b>Calibrator replicates</b>	<b>3</b>
	<b>Blank replicates</b>	<b>3</b>
	<b>Calibration curve</b>	<b>-</b>
<b>OPTIONS</b>	<b>Blank absorbance limit</b>	<b>0.700</b>
	<b>Kinetic blank limit</b>	<b>-</b>
	<b>Linerity limit</b>	<b>800</b>

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<b>GENERAL</b>	<b>Test name</b>	<b>Albumina</b>
	<b>Analysis mode</b>	<b>Endpoint mon.</b>
	<b>Sample type</b>	<b>Serum</b>
	<b>Units</b>	<b>g/dL</b>
	<b>Reaction type</b>	<b>Increasing</b>
	<b>Decimals</b>	<b>2</b>
	<b>Replicates</b>	<b>1</b>
	<b>Name of assoc. constituent</b>	<b>-</b>
<b>PROCEDURE</b>	<b>Type of reading</b>	<b>monochrom.</b>
	<b>Volumes</b>	
	<b>Sample</b>	<b>3</b>
	<b>Reagent 1</b>	<b>300</b>
	<b>Reagent 2</b>	<b>-</b>
	<b>Washing</b>	<b>1.2</b>
	<b>Predilution factor</b>	<b>-</b>
	<b>Filters</b>	
	<b>Main</b>	<b>635</b>
	<b>Reference</b>	
	<b>Times</b>	
	<b>Reading 1</b>	<b>120 s</b>
<b>Reading 2</b>	<b>-</b>	
<b>Reagent 2</b>	<b>-</b>	
<b>Postdilution factor</b>	<b>2</b>	
<b>CALIBRATION</b>	<b>Type of calibration</b>	<b>Múltiplo (*)</b>
	<b>Calibrator replicates</b>	<b>3</b>
	<b>Blank replicates</b>	<b>3</b>
	<b>Calibration curve</b>	<b>-</b>
<b>OPTIONS</b>	<b>Blank absorbance limit</b>	<b>0.200</b>
	<b>Kinetic blank limit</b>	<b>-</b>
	<b>Linerity limit</b>	<b>6.0</b>

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GENERAL	Test name	ALT
	Analysis mode	Kinetic mon.
	Sample type	Serum
	Units	U/L
	Reaction type	decreasing
	Decimals	0
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monoch.
	Volumes	
	Sample	15
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	340
	Reference	-
	Times	
Reading 1	96 s	
Reading 2	264	
Reagent 2	-	
Postdilution factor	2	
CALIBRATION	Type of calibration	Múltiplo (*)
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	1.100
	Kinetic blank limit	-
	Linerity limit	400

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<b>GENERAL</b>	<b>Test name</b>	<b>Amilase</b>
	<b>Analysis mode</b>	<b>Kinetic mon.</b>
	<b>Sample type</b>	<b>Serum</b>
	<b>Units</b>	<b>U/L</b>
	<b>Reaction type</b>	<b>Increasing</b>
	<b>Decimals</b>	<b>0</b>
	<b>Replicates</b>	<b>1</b>
	<b>Name of assoc. constituent</b>	<b>-</b>
<b>PROCEDURE</b>	<b>Type of reading</b>	<b>Monoch.</b>
	<b>Volumes</b>	
	<b>Sample</b>	<b>6</b>
	<b>Reagent 1</b>	<b>300</b>
	<b>Reagent 2</b>	<b>-</b>
	<b>Washing</b>	<b>1.2</b>
	<b>Predilution factor</b>	<b>-</b>
	<b>Filters</b>	
	<b>Main</b>	<b>405</b>
	<b>Reference</b>	<b>-</b>
<b>Times</b>	<b>Reading 1</b>	<b>72 s</b>
	<b>Reading 2</b>	<b>216 s</b>
	<b>Reagent 2</b>	<b>-</b>
<b>CALIBRATION</b>	<b>Postdilution factor</b>	<b>2</b>
	<b>Type of calibration</b>	<b>Múltiplo (*)</b>
	<b>Calibrator replicates</b>	<b>3</b>
	<b>Blank replicates</b>	<b>3</b>
	<b>Calibration curve</b>	<b>-</b>
<b>OPTIONS</b>	<b>Blank absorbance limit</b>	<b>0.700</b>
	<b>Kinetic blank limit</b>	<b>-</b>
	<b>Linerity limit</b>	<b>2000</b>

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GENERAL	Test name	AST	
	Analysis mode	kinetic mon.	
	Sample type	Serum	
	Units	U/L	
	Reaction type	decreasing	
	Decimals	0	
	Replicates	1	
Name of assoc. constituent		-	
PROCEDURE	Volumes	Type of reading	Monoch.
		Sample	15
		Reagent 1	300
		Reagent 2	-
		Washing	1.2
	Filters	Predilution factor	-
		Main	340
		Reference	-
	Times	Reading 1	96 s
		Reading 2	264 s
		Reagent 2	-
Postdilution factor		2	
CALIBRATION	Type of calibration	Múltiplo (*)	
	Calibrator replicates	3	
	Blank replicates	3	
	Calibration curve	-	
OPTIONS	Blank absorbance limit	1.100	
	Kinetic blank limit	-	
	Linerity limit	400	

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<b>GENERAL</b>	<b>Test name</b>	<b>Cálcio CPC</b>	
	<b>Analysis mode</b>	<b>Endpoint mon.</b>	
	<b>Sample type</b>	<b>Serum</b>	
	<b>Units</b>	<b>mg/dL</b>	
	<b>Reaction type</b>	<b>Increasing</b>	
	<b>Decimals</b>	<b>1</b>	
	<b>Replicates</b>	<b>1</b>	
	<b>Name of assoc. constituent</b>	<b>-</b>	
<b>PROCEDURE</b>	<b>Volumes</b>	<b>Type of reading</b>	<b>Monoch.</b>
		<b>Sample</b>	<b>3</b>
		<b>Reagent 1</b>	<b>300</b>
		<b>Reagent 2</b>	<b>----</b>
		<b>Washing</b>	<b>1.2</b>
		<b>Predilution factor</b>	<b>-</b>
	<b>Filters</b>	<b>Main</b>	<b>560</b>
		<b>Reference</b>	<b>-</b>
	<b>Times</b>	<b>Reading 1</b>	<b>30 s</b>
		<b>Reading 2</b>	<b>----</b>
		<b>Reagent 2</b>	<b>----</b>
		<b>Postdilution factor</b>	<b>2</b>
	<b>CALIBRATION</b>	<b>Type of calibration</b>	<b>Múltiplo (*)</b>
<b>Calibrator replicates</b>		<b>3</b>	
<b>Blank replicates</b>		<b>3</b>	
<b>Calibration curve</b>		<b>-</b>	
<b>OPTIONS</b>	<b>Blank absorbance limit</b>	<b>2.000</b>	
	<b>Kinetic blank limit</b>	<b>-</b>	
	<b>Linerity limit</b>	<b>16.0</b>	

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GENERAL	Test name	Cálcio Arsenaso
	Analysis mode	Endpoint mon.
	Sample type	Serum
	Units	mg/dL
	Reaction type	Increasing
	Decimals	1
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monoch.
	Volumes	
	Sample	3
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	670
	Reference	-
	Times	
	Reading 1	120 s
Reading 2	-	
Reagent 2	-	
Postdilution factor	2	
CALIBRATION	Type of calibration	Múltiplo (*)
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	1.000
	Kinetic blank limit	-
	Linerity limit	17

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<b>GENERAL</b>	<b>Test name</b>	<b>CLF</b>	
	<b>Analysis mode</b>	<b>Differential Bir.</b>	
	<b>Sample type</b>	<b>Serum</b>	
	<b>Units</b>	<b>µg/dL</b>	
	<b>Reaction type</b>	<b>decreasing</b>	
	<b>Decimals</b>	<b>1</b>	
	<b>Replicates</b>	<b>1</b>	
	<b>Name of assoc. constituent</b>	<b>-</b>	
<b>PROCEDURE</b>	<b>Volumes</b>	<b>Type of reading</b>	<b>Monoch.</b>
		<b>Sample</b>	<b>30</b>
		<b>Reagent 1</b>	<b>235</b>
		<b>Reagent 2</b>	<b>65</b>
		<b>Washing</b>	<b>1.2</b>
		<b>Predilution factor</b>	<b>-</b>
	<b>Filters</b>	<b>Main</b>	<b>560</b>
		<b>Reference</b>	<b>-</b>
		<b>Times</b>	<b>Reading 1</b>
	<b>Reading 2</b>		<b>600 s</b>
	<b>Reagent 2</b>		<b>312 s</b>
	<b>Postdilution factor</b>		<b>2</b>
	<b>CALIBRATION</b>	<b>Type of calibration</b>	<b>specific (*)</b>
<b>Calibrator replicates</b>		<b>2</b>	
<b>Blank replicates</b>		<b>2</b>	
<b>Calibration curve</b>		<b>-</b>	
<b>OPTIONS</b>		<b>Blank absorbance limit</b>	<b>0.350</b>
	<b>Kinetic blank limit</b>	<b>-</b>	
	<b>Linerity limit</b>	<b>500</b>	

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GENERAL	Test name	CKMB
	Analysis mode	Kinetic mon.
	Sample type	Serum
	Units	U/L
	Reaction type	Increasing
	Decimals	0
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monoch.
	Volumes	
	Sample	15
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	340
	Reference	-
	Times	
	Reading 1	312 s
Reading 2	600 s	
Reagent 2	-	
Postdilution factor	-	
CALIBRATION	Type of calibration	Specific (*)
	Factor	
	Calibrator replicates	
	Blank replicates	3
OPTIONS	Calibration curve	-
	Blank absorbance limit	0.500
	Kinetic blank limit	-
	Linerity limit	500

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GENERAL	Test name	CK NAC
	Analysis mode	Kinetic mon.
	Sample type	Serum
	Units	U/L
	Reaction type	Increasing
	Decimals	0
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monoch.
	Volumes	
	Sample	6
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	340
	Reference	-
	Times	
	Reading 1	120 s
Reading 2	240 s	
Reagent 2	-	
Postdilution factor	2	
CALIBRATION	Type of calibration	specific (*)
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	0.500
	Kinetic blank limit	-
	Linerity limit	2000

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<b>GENERAL</b>	<b>Test name</b>	<b>Colesterol</b>	
	<b>Analysis mode</b>	<b>Endpoint mon.</b>	
	<b>Sample type</b>	<b>Serum</b>	
	<b>Units</b>	<b>mg/dL</b>	
	<b>Reaction type</b>	<b>Increasing</b>	
	<b>Decimals</b>	<b>0</b>	
	<b>Replicates</b>	<b>1</b>	
	<b>Name of assoc. constituent</b>	<b>-</b>	
<b>PROCEDURE</b>	<b>Volumes</b>	<b>Type of reading</b>	<b>Bichrom.</b>
		<b>Sample</b>	<b>3 µL</b>
		<b>Reagent 1</b>	<b>300 µL</b>
		<b>Reagent 2</b>	<b>-</b>
		<b>Washing</b>	<b>1.2</b>
	<b>Filters</b>	<b>Predilution factor</b>	<b>-</b>
		<b>Main</b>	<b>505</b>
	<b>Times</b>	<b>Reference</b>	<b>670</b>
		<b>Reading 1</b>	<b>600 s</b>
		<b>Reading 2</b>	<b>-</b>
		<b>Reagent 2</b>	<b>-</b>
		<b>Postdilution factor</b>	<b>-</b>
<b>CALIBRATION</b>	<b>Type of calibration</b>	<b>Múltiplo (*)</b>	
	<b>Calibrator replicates</b>	<b>1</b>	
	<b>Blank replicates</b>	<b>1</b>	
	<b>Calibration curve</b>	<b>-</b>	
<b>OPTIONS</b>	<b>Blank absorbance limit</b>	<b>0.180</b>	
	<b>Kinetic blank limit</b>	<b>-</b>	
	<b>Linerity limit</b>	<b>500 mg/dL</b>	

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<b>GENERAL</b>	<b>Test name</b>	<b>Colesterol HDL</b>	
	<b>Analysis mode</b>	<b>Endpoint mon.</b>	
	<b>Sample type</b>	<b>Serum</b>	
	<b>Units</b>	<b>mg/dL</b>	
	<b>Reaction type</b>	<b>Increasing</b>	
	<b>Decimals</b>	<b>0</b>	
	<b>Replicates</b>	<b>1</b>	
	<b>Name of assoc. constituent</b>	<b>-</b>	
<b>PROCEDURE</b>	<b>Volumes</b>	<b>Type of reading</b>	<b>Bichrom.</b>
		<b>Sample</b>	<b>30 µL</b>
		<b>Reagent 1</b>	<b>300 µL</b>
		<b>Reagent 2</b>	<b>-</b>
		<b>Washing</b>	<b>1.2</b>
		<b>Predilution factor</b>	<b>-</b>
	<b>Filters</b>	<b>Main</b>	<b>505</b>
		<b>Reference</b>	<b>670</b>
	<b>Times</b>	<b>Reading 1</b>	<b>600 s</b>
		<b>Reading 2</b>	<b>-</b>
		<b>Reagent 2</b>	<b>-</b>
		<b>Postdilution factor</b>	<b>2</b>
<b>CALIBRATION</b>	<b>Type of calibration</b>	<b>specific (*)</b>	
	<b>Calibrator replicates</b>	<b>1</b>	
	<b>Blank replicates</b>	<b>1</b>	
	<b>Calibration curve</b>	<b>-</b>	
<b>OPTIONS</b>	<b>Blank absorbance limit</b>	<b>0.180</b>	
	<b>Kinetic blank limit</b>	<b>-</b>	
	<b>Linerity limit</b>	<b>200 mg/dL</b>	

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<b>GENERAL</b>	<b>Test name</b>	<b>CHE</b>
	<b>Analysis mode</b>	<b>Kinetic bir.</b>
	<b>Sample type</b>	<b>Serum</b>
	<b>Units</b>	<b>U/L</b>
	<b>Reaction type</b>	<b>Decreasing</b>
	<b>Decimals</b>	<b>0</b>
	<b>Replicates</b>	<b>1</b>
	<b>Name of assoc. constituent</b>	<b>-</b>
<b>PROCEDURE</b>	<b>Type of reading</b>	<b>Monoch.</b>
	<b>Volumes</b>	
	<b>Sample</b>	<b>4</b>
	<b>Reagent 1</b>	<b>200</b>
	<b>Reagent 2</b>	<b>50</b>
	<b>Washing</b>	<b>1.2</b>
	<b>Predilution factor</b>	<b>-</b>
	<b>Filters</b>	
	<b>Main</b>	<b>405</b>
	<b>Reference</b>	<b>-</b>
	<b>Times</b>	
	<b>Reading 1</b>	<b>312 s</b>
<b>Reading 2</b>	<b>408 s</b>	
<b>Reagent 2</b>	<b>192 s</b>	
<b>Postdilution factor</b>	<b>2</b>	
<b>CALIBRATION</b>	<b>Type of calibration</b>	<b>Multiple (*)</b>
	<b>Factor</b>	
	<b>Calibrator replicates</b>	<b>3</b>
	<b>Blank replicates</b>	<b>3</b>
<b>OPTIONS</b>	<b>Calibration curve</b>	<b>-</b>
	<b>Blank absorbance limit</b>	<b>1.100</b>
	<b>Kinetic blank limit</b>	<b>-</b>
	<b>Linerity limit</b>	<b>20000</b>

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<b>GENERAL</b>	<b>Test name</b>	<b>Creatinina</b>	
	<b>Analysis mode</b>	<b>Fixed-time mon.</b>	
	<b>Sample type</b>	<b>Serum</b>	
	<b>Units</b>	<b>mg/dL</b>	
	<b>Reaction type</b>	<b>Increasing</b>	
	<b>Decimals</b>	<b>2</b>	
	<b>Replicates</b>	<b>1</b>	
	<b>Name of assoc. constituent</b>	<b>-</b>	
<b>PROCEDURE</b>	<b>Volumes</b>	<b>Type of reading</b>	<b>Monoch.</b>
		<b>Sample</b>	<b>30 µL</b>
		<b>Reagent 1</b>	<b>300 µL</b>
		<b>Reagent 2</b>	<b>-</b>
		<b>Washing</b>	<b>1.2</b>
	<b>Filters</b>	<b>Predilution factor</b>	<b>-</b>
		<b>Main</b>	<b>505</b>
		<b>Reference</b>	<b>-</b>
	<b>Times</b>	<b>Reading 1</b>	<b>48 s</b>
		<b>Reading 2</b>	<b>120 s</b>
		<b>Reagent 2</b>	<b>-</b>
		<b>Postdilution factor</b>	<b>2</b>
<b>CALIBRATION</b>	<b>Type of calibration</b>	<b>Múltiplo (*)</b>	
	<b>Calibrator replicates</b>	<b>3</b>	
	<b>Blank replicates</b>	<b>3</b>	
	<b>Calibration curve</b>	<b>-</b>	
<b>OPTIONS</b>	<b>Blank absorbance limit</b>	<b>0.200</b>	
	<b>Kinetic blank limit</b>	<b>-</b>	
	<b>Linerity limit</b>	<b>12.0 mg/dL</b>	

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<b>GENERAL</b>	<b>Test name</b>	<b>Creatinina</b>
	<b>Analysis mode</b>	<b>Fixed-time mon.</b>
	<b>Sample type</b>	<b>Serum</b>
	<b>Units</b>	<b>mg/dL</b>
	<b>Reaction type</b>	<b>Increasing</b>
	<b>Decimals</b>	<b>2</b>
	<b>Replicates</b>	<b>1</b>
	<b>Name of assoc. constituent</b>	<b>-</b>
<b>PROCEDURE</b>	<b>Type of reading</b>	<b>Monoch.</b>
	<b>Volumes</b>	
	<b>Sample</b>	<b>15 µL</b>
	<b>Reagent 1</b>	<b>300 µL</b>
	<b>Reagent 2</b>	<b>-</b>
	<b>Washing</b>	<b>1.2</b>
	<b>Predilution factor</b>	<b>-</b>
	<b>Filters</b>	
	<b>Main</b>	<b>505</b>
	<b>Reference</b>	<b>-</b>
	<b>Times</b>	
	<b>Reading 1</b>	<b>72 s</b>
<b>Reading 2</b>	<b>240 s</b>	
	<b>Reagent 2</b>	<b>-</b>
	<b>Postdilution factor</b>	<b>2</b>
<b>CALIBRATION</b>	<b>Type of calibration</b>	<b>Múltiplo (*)</b>
	<b>Calibrator replicates</b>	<b>2</b>
	<b>Blank replicates</b>	<b>1</b>
	<b>Calibration curve</b>	<b>-</b>
<b>OPTIONS</b>	<b>Blank absorbance limit</b>	<b>0.200</b>
	<b>Kinetic blank limit</b>	<b>-</b>
	<b>Linerity limit</b>	<b>12 mg/dL</b>

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<b>GENERAL</b>	<b>Test name</b>	<b>Ferritina</b>	
	<b>Analysis mode</b>	<b>Fixed Time mon.</b>	
	<b>Sample type</b>	<b>Serum</b>	
	<b>Units</b>	<b>µg/L</b>	
	<b>Reaction type</b>	<b>Increasing</b>	
	<b>Turbidim. Test</b>	<b>yes</b>	
	<b>Decimals</b>	<b>0</b>	
	<b>Replicates</b>	<b>1</b>	
	<b>Name of assoc. constituent</b>	<b>-</b>	
<b>PROCEDURE</b>	<b>Type of reading</b>	<b>monoch.</b>	
	<b>Volumes</b>	<b>Sample</b>	<b>8</b>
		<b>Reagent 1</b>	<b>240</b>
		<b>Reagent 2</b>	<b>-</b>
		<b>Washing</b>	<b>1.2</b>
		<b>Predilution factor</b>	<b>-</b>
	<b>Filters</b>	<b>Main</b>	<b>535</b>
		<b>Reference</b>	<b>-</b>
	<b>Times</b>	<b>Reading 1</b>	<b>48 s</b>
		<b>Reading 2</b>	<b>360 s</b>
		<b>Reagent 2</b>	<b>-</b>
		<b>Postdilution factor</b>	<b>-</b>
<b>CALIBRATION</b>	<b>Type of calibration</b>	<b>specific (*)</b>	
	<b>Calibrator number</b>	<b>5</b>	
	<b>Calibrator replicates</b>	<b>3</b>	
	<b>Blank replicates</b>	<b>3</b>	
	<b>Calibration curve</b>	<b>Increasing Polygonal</b>	
<b>OPTIONS</b>	<b>Blank absorbance limit</b>	<b>1.600</b>	
	<b>Kinetic blank limit</b>	<b>-</b>	
	<b>Linerity limit</b>	<b>-</b>	

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<b>GENERAL</b>	<b>Test name</b>	<b>Ferro Ferrozina</b>
	<b>Analysis mode</b>	<b>Differential bir.</b>
	<b>Sample type</b>	<b>Serum</b>
	<b>Units</b>	<b>µg/dL</b>
	<b>Reaction type</b>	<b>Increasing</b>
	<b>Decimals</b>	<b>0</b>
	<b>Replicates</b>	<b>1</b>
	<b>Name of assoc. constituent</b>	<b>-</b>
<b>PROCEDURE</b>	<b>Type of reading</b>	<b>Monoch.</b>
	<b>Volumes</b>	
	<b>Sample</b>	<b>25 µL</b>
	<b>Reagent 1</b>	<b>200 µL</b>
	<b>Reagent 2</b>	<b>50 µL</b>
	<b>Washing</b>	<b>1.2</b>
	<b>Predilution factor</b>	<b>-</b>
	<b>Filters</b>	
	<b>Main</b>	<b>560</b>
	<b>Reference</b>	<b>-</b>
	<b>Times</b>	
	<b>Reading 1</b>	<b>288 s</b>
<b>Reading 2</b>	<b>600 s</b>	
<b>Reagent 2</b>	<b>312 s</b>	
<b>Postdilution factor</b>	<b>2</b>	
<b>CALIBRATION</b>	<b>Type of calibration</b>	<b>Específico (*)</b>
	<b>Calibrator replicates</b>	<b>2</b>
	<b>Blank replicates</b>	<b>2</b>
	<b>Calibration curve</b>	<b>-</b>
<b>OPTIONS</b>	<b>Blank absorbance limit</b>	<b>0.300</b>
	<b>Kinetic blank limit</b>	<b>-</b>
	<b>Linerity limit</b>	<b>1000 µg/dL</b>

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<b>GENERAL</b>	<b>Test name</b>	<b>Ferro Cromazurol</b>	
	<b>Analysis mode</b>	<b>Endpoint mon.</b>	
	<b>Sample type</b>	<b>Serum</b>	
	<b>Units</b>	<b>µg/dL</b>	
	<b>Reaction type</b>	<b>Increasing</b>	
	<b>Decimals</b>	<b>0</b>	
	<b>Replicates</b>	<b>1</b>	
	<b>Name of assoc. constituent</b>	<b>-</b>	
<b>PROCEDURE</b>	<b>Volumes</b>	<b>Type of reading</b>	<b>Monoch.</b>
		<b>Sample</b>	<b>15</b>
		<b>Reagent 1</b>	<b>300</b>
		<b>Reagent 2</b>	<b>-</b>
		<b>Washing</b>	<b>1.2</b>
		<b>Predilution factor</b>	<b>-</b>
	<b>Filters</b>	<b>Main</b>	<b>635</b>
		<b>Reference</b>	<b>-</b>
	<b>Times</b>	<b>Reading 1</b>	<b>312 s</b>
		<b>Reading 2</b>	<b>-</b>
		<b>Reagent 2</b>	<b>-</b>
		<b>Postdilution factor</b>	<b>2</b>
<b>CALIBRATION</b>	<b>Type of calibration</b>	<b>Múltiplo (*)</b>	
	<b>Calibrator replicates</b>	<b>3</b>	
	<b>Blank replicates</b>	<b>3</b>	
	<b>Calibration curve</b>	<b>-</b>	
<b>OPTIONS</b>	<b>Blank absorbance limit</b>	<b>0.750</b>	
	<b>Kinetic blank limit</b>	<b>-</b>	
	<b>Linerity limit</b>	<b>500</b>	

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<b>GENERAL</b>	<b>Test name</b>	<b>Fosfatase Alcalina AMP</b>	
	<b>Analysis mode</b>	<b>Kinetic mon.</b>	
	<b>Sample type</b>	<b>Serum</b>	
	<b>Units</b>	<b>U/L</b>	
	<b>Reaction type</b>	<b>Increasing</b>	
	<b>Decimals</b>	<b>0</b>	
	<b>Replicates</b>	<b>1</b>	
	<b>Name of assoc. constituent</b>	<b>-</b>	
<b>PROCEDURE</b>	<b>Volumes</b>	<b>Type of reading</b>	<b>Monoch.</b>
		<b>Sample</b>	<b>6</b>
		<b>Reagent 1</b>	<b>300</b>
		<b>Reagent 2</b>	<b>-</b>
		<b>Washing</b>	<b>1.2</b>
		<b>Predilution factor</b>	<b>-</b>
	<b>Filters</b>	<b>Main</b>	<b>405</b>
		<b>Reference</b>	<b>-</b>
	<b>Times</b>	<b>Reading 1</b>	<b>72 s</b>
		<b>Reading 2</b>	<b>192 s</b>
		<b>Reagent 2</b>	<b>-</b>
		<b>Postdilution factor</b>	<b>2</b>
<b>CALIBRATION</b>	<b>Type of calibration</b>	<b>Múltiplo (*)</b>	
	<b>Calibrator replicates</b>	<b>3</b>	
	<b>Blank replicates</b>	<b>3</b>	
	<b>Calibration curve</b>	<b>-</b>	
<b>OPTIONS</b>	<b>Blank absorbance limit</b>	<b>1.000</b>	
	<b>Kinetic blank limit</b>	<b>-</b>	
	<b>Linerity limit</b>	<b>1500</b>	

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<b>GENERAL</b>		<b>Test name</b>	<b>Fósforo UV</b>
		<b>Analysis mode</b>	<b>Endpoint mon.</b>
		<b>Sample type</b>	<b>Serum</b>
		<b>Units</b>	<b>mg/dL</b>
		<b>Reaction type</b>	<b>Increasing</b>
		<b>Decimals</b>	<b>2</b>
		<b>Replicates</b>	<b>1</b>
		<b>Name of assoc. constituent</b>	<b>-</b>
<b>PROCEDURE</b>	<b>Volumes</b>	<b>Type of reading</b>	<b>Monoch.</b>
		<b>Sample</b>	<b>3</b>
		<b>Reagent 1</b>	<b>300</b>
		<b>Reagent 2</b>	<b>-</b>
		<b>Washing</b>	<b>1.2</b>
	<b>Filters</b>	<b>Predilution factor</b>	<b>-</b>
		<b>Main</b>	<b>340</b>
	<b>Times</b>	<b>Reference</b>	<b>-</b>
		<b>Reading 1</b>	<b>312 s</b>
		<b>Reading 2</b>	<b>-</b>
		<b>Reagent 2</b>	<b>-</b>
		<b>Postdilution factor</b>	<b>2</b>
<b>CALIBRATION</b>	<b>Type of calibration</b>	<b>Múltiplo (*)</b>	
	<b>Calibrator replicates</b>	<b>1</b>	
	<b>Blank replicates</b>	<b>1</b>	
	<b>Calibration curve</b>	<b>-</b>	
<b>OPTIONS</b>	<b>Blank absorbance limit</b>	<b>0.500</b>	
	<b>Kinetic blank limit</b>	<b>-</b>	
	<b>Linerity limit</b>	<b>20.0 mg/dL</b>	

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<b>GENERAL</b>	<b>Test name</b>	<b>FR Turbidim</b>
	<b>Analysis mode</b>	<b>Endpoint bir.</b>
	<b>Sample type</b>	<b>Serum</b>
	<b>Units</b>	<b>UI/mL</b>
	<b>Reaction type</b>	<b>Increasing</b>
	<b>Turbidim. Test</b>	<b>yes</b>
	<b>Decimals</b>	<b>0</b>
	<b>Replicates</b>	<b>1</b>
	<b>Name of assoc. constituent</b>	<b>-</b>
<b>PROCEDURE</b>	<b>Type of reading</b>	<b>monoch.</b>
	<b>Volumes</b>	
	<b>Sample</b>	<b>3</b>
	<b>Reagent 1</b>	<b>240</b>
	<b>Reagent 2</b>	<b>60</b>
	<b>Washing</b>	<b>1.2</b>
	<b>Predilution factor</b>	<b>-</b>
	<b>Filters</b>	
	<b>Main</b>	<b>635</b>
	<b>Reference</b>	<b>-</b>
	<b>Times</b>	
	<b>Reading 1</b>	<b>144 s</b>
<b>Reading 2</b>	<b>-</b>	
<b>Reagent 2</b>	<b>24 s</b>	
<b>Postdilution factor</b>	<b>2</b>	
<b>CALIBRATION</b>	<b>Type of calibration</b>	<b>specific (*)</b>
	<b>Calibrator number</b>	<b>5</b>
	<b>Calibrator replicates</b>	<b>3</b>
	<b>Blank replicates</b>	<b>3</b>
	<b>Calibration curve</b>	<b>Increasing Polygonal</b>
<b>OPTIONS</b>	<b>Blank absorbance limit</b>	<b>1.400</b>
	<b>Kinetic blank limit</b>	<b>-</b>
	<b>Linerity limit</b>	<b>-</b>

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<b>GENERAL</b>	<b>Test name</b>	<b>Frutosamina</b>	
	<b>Analysis mode</b>	<b>Fixed-time mon.</b>	
	<b>Sample type</b>	<b>Serum</b>	
	<b>Units</b>	<b>mmol/L</b>	
	<b>Reaction type</b>	<b>Increasing</b>	
	<b>Decimals</b>	<b>2</b>	
	<b>Replicates</b>	<b>1</b>	
	<b>Name of assoc. constituent</b>	<b>-</b>	
<b>PROCEDURE</b>	<b>Volumes</b>	<b>Type of reading</b>	<b>Monoch.</b>
		<b>Sample</b>	<b>15</b>
		<b>Reagent 1</b>	<b>300</b>
		<b>Reagent 2</b>	<b>-</b>
		<b>Washing</b>	<b>1.2</b>
		<b>Predilution factor</b>	<b>-</b>
	<b>Filters</b>	<b>Main</b>	<b>535</b>
		<b>Reference</b>	<b>-</b>
	<b>Times</b>	<b>Reading 1</b>	<b>408 s</b>
		<b>Reading 2</b>	<b>600 s</b>
		<b>Reagent 2</b>	<b>-</b>
		<b>Postdilution factor</b>	<b>2</b>
	<b>CALIBRATION</b>	<b>Type of calibration</b>	<b>Specific (*)</b>
<b>Calibrator replicates</b>		<b>3</b>	
<b>Blank replicates</b>		<b>3</b>	
<b>Calibration curve</b>		<b>-</b>	
<b>Blank absorbance limit</b>		<b>0.065</b>	
<b>OPTIONS</b>	<b>Kinetic blank limit</b>	<b>-</b>	
	<b>Linerity limit</b>	<b>7.0</b>	

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GENERAL	Test name	Gama GT
	Analysis mode	Kinetic mon.
	Sample type	Serum
	Units	U/L
	Reaction type	Increasing
	Decimals	0
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monoch.
	Sample	15
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Main	405
	Reference	-
	Reading 1	72 s
	Reading 2	192 s
	Reagent 2	-
	Postdilution factor	2
CALIBRATION	Type of calibration	Múltiplo (*)
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	1.450
	Kinetic blank limit	-
	Linerity limit	700

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<b>GENERAL</b>	<b>Test name</b>	<b>Glicose</b>
	<b>Analysis mode</b>	<b>Endpoint mon.</b>
	<b>Sample type</b>	<b>Serum</b>
	<b>Units</b>	<b>mg/dL</b>
	<b>Reaction type</b>	<b>Increasing</b>
	<b>Decimals</b>	<b>0</b>
	<b>Replicates</b>	<b>1</b>
	<b>Name of assoc. constituent</b>	<b>-</b>
<b>PROCEDURE</b>	<b>Type of reading</b>	<b>bichrom.</b>
	<b>Volumes</b>	
	<b>Sample</b>	<b>3</b>
	<b>Reagent 1</b>	<b>300</b>
	<b>Reagent 2</b>	<b>-</b>
	<b>Washing</b>	<b>1.2</b>
	<b>Predilution factor</b>	<b>-</b>
	<b>Filters</b>	
	<b>Main</b>	<b>505</b>
	<b>Reference</b>	<b>670</b>
	<b>Times</b>	
	<b>Reading 1</b>	<b>312 s</b>
<b>Reading 2</b>	<b>-</b>	
<b>Reagent 2</b>	<b>-</b>	
<b>Postdilution factor</b>	<b>2</b>	
<b>CALIBRATION</b>	<b>Type of calibration</b>	<b>Múltiplo (*)</b>
	<b>Calibrator replicates</b>	<b>1</b>
	<b>Blank replicates</b>	<b>1</b>
	<b>Calibration curve</b>	<b>-</b>
<b>OPTIONS</b>	<b>Blank absorbance limit</b>	<b>0.180</b>
	<b>Kinetic blank limit</b>	<b>-</b>
	<b>Linerity limit</b>	<b>500 mg/dL</b>

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<b>GENERAL</b>	<b>Test name</b>	<b>HDL Direto</b>
	<b>Analysis mode</b>	<b>Diferencial birreagente</b>
	<b>Sample type</b>	<b>Serum</b>
	<b>Units</b>	<b>mg/dL</b>
	<b>Reaction type</b>	<b>Increasing</b>
	<b>Decimals</b>	<b>0</b>
	<b>Replicates</b>	<b>1</b>
	<b>Name of assoc. constituent</b>	<b>-</b>
<b>PROCEDURE</b>	<b>Type of reading</b>	<b>Monocromática</b>
	<b>Volumes</b>	
	<b>Sample</b>	<b>3</b>
	<b>Reagent 1</b>	<b>300</b>
	<b>Reagent 2</b>	<b>100</b>
	<b>Washing</b>	<b>1.2</b>
	<b>Predilution factor</b>	<b>-</b>
	<b>Filters</b>	
	<b>Main</b>	<b>535</b>
	<b>Reference</b>	
	<b>Times</b>	
	<b>Reading 1</b>	<b>168 s</b>
<b>Reading 2</b>	<b>480 s</b>	
<b>Reagent 2</b>	<b>192 s</b>	
	<b>Postdilution factor</b>	<b>2</b>
<b>CALIBRATION</b>	<b>Type of calibration</b>	<b>Specific (*)</b>
	<b>Calibrator replicates</b>	<b>1</b>
	<b>Blank replicates</b>	<b>2</b>
	<b>Calibration curve</b>	<b>-</b>
<b>OPTIONS</b>	<b>Blank absorbance limit</b>	<b>0.200</b>
	<b>Kinetic blank limit</b>	<b>-</b>
	<b>Linerity limit</b>	<b>200 mg/dL</b>

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<b>GENERAL</b>	<b>Test name</b>	<b>LDL Direto</b>
	<b>Analysis mode</b>	<b>Differential bir.</b>
	<b>Sample type</b>	<b>Serum</b>
	<b>Units</b>	<b>mg/dL</b>
	<b>Reaction type</b>	<b>Increasing</b>
	<b>Decimals</b>	<b>0</b>
	<b>Replicates</b>	<b>1</b>
	<b>Name of assoc. constituent</b>	<b>-</b>
<b>PROCEDURE</b>	<b>Type of reading</b>	<b>Monoch.</b>
	<b>Volumes</b>	
	<b>Sample</b>	<b>3</b>
	<b>Reagent 1</b>	<b>300</b>
	<b>Reagent 2</b>	<b>100</b>
	<b>Washing</b>	<b>1.2</b>
	<b>Predilution factor</b>	<b>-</b>
	<b>Filters</b>	
	<b>Main</b>	<b>535</b>
	<b>Reference</b>	<b>-</b>
	<b>Times</b>	
	<b>Reading 1</b>	<b>288 s</b>
<b>Reading 2</b>	<b>480 s</b>	
<b>Reagent 2</b>	<b>312 s</b>	
<b>Postdilution factor</b>	<b>2</b>	
<b>CALIBRATION</b>	<b>Type of calibration</b>	<b>specific (*)</b>
	<b>Calibrator replicates</b>	<b>3</b>
	<b>Blank replicates</b>	<b>3</b>
	<b>Calibration curve</b>	<b>-</b>
<b>OPTIONS</b>	<b>Blank absorbance limit</b>	<b>0.100</b>
	<b>Kinetic blank limit</b>	<b>-</b>
	<b>Linerity limit</b>	<b>990</b>

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GENERAL	Test name	LDH
	Analysis mode	Kinetic mon.
	Sample type	Serum
	Units	U/L
	Reaction type	decreasing
	Decimals	0
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	Monoch.
	Sample	6
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Main	340
	Reference	-
	Reading 1	72 s
	Reading 2	192 s
Reagent 2	-	
Postdilution factor	2	
CALIBRATION	Type of calibration	Múltiplo (*)
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	0.800
	Kinetic blank limit	-
	Linerity limit	1500

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<b>GENERAL</b>	<b>Test name</b>	<b>LIPASE</b>
	<b>Analysis mode</b>	<b>Fixed time Bir.</b>
	<b>Sample type</b>	<b>Serum</b>
	<b>Units</b>	<b>µg/dL</b>
	<b>Reaction type</b>	<b>Increasing</b>
	<b>Decimals</b>	<b>1</b>
	<b>Replicates</b>	<b>1</b>
	<b>Name of assoc. constituent</b>	<b>-</b>
<b>PROCEDURE</b>	<b>Type of reading</b>	<b>Monoch.</b>
	<b>Volumes</b>	
	<b>Sample</b>	<b>3</b>
	<b>Reagent 1</b>	<b>200</b>
	<b>Reagent 2</b>	<b>110</b>
	<b>Washing</b>	<b>1.2</b>
	<b>Predilution factor</b>	<b>-</b>
	<b>Filters</b>	
	<b>Main</b>	<b>560</b>
	<b>Reference</b>	<b>-</b>
	<b>Times</b>	
	<b>Reading 1</b>	<b>120 s</b>
	<b>Reading 2</b>	<b>216 s</b>
<b>Reagent 2</b>	<b>24 s</b>	
<b>Postdilution factor</b>	<b>2</b>	
<b>CALIBRATION</b>	<b>Type of calibration</b>	<b>multiple (*)</b>
	<b>Calibrator replicates</b>	<b>2</b>
	<b>Blank replicates</b>	<b>2</b>
	<b>Calibration curve</b>	<b>-</b>
<b>OPTIONS</b>	<b>Blank absorbance limit</b>	<b>0.200</b>
	<b>Kinetic blank limit</b>	<b>-</b>
	<b>Linerity limit</b>	<b>300</b>

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<b>GENERAL</b>	<b>Test name</b>	<b>Magnésio</b>
	<b>Analysis mode</b>	<b>Endpoint mon.</b>
	<b>Sample type</b>	<b>Serum</b>
	<b>Units</b>	<b>mg/dL</b>
	<b>Reaction type</b>	<b>Increasing</b>
	<b>Decimals</b>	<b>2</b>
	<b>Replicates</b>	<b>1</b>
	<b>Name of assoc. constituent</b>	<b>-</b>
<b>PROCEDURE</b>	<b>Type of reading</b>	<b>Monoch.</b>
	<b>Volumes</b>	
	<b>Sample</b>	<b>3</b>
	<b>Reagent 1</b>	<b>300</b>
	<b>Reagent 2</b>	<b>-</b>
	<b>Washing</b>	<b>1.2</b>
	<b>Predilution factor</b>	<b>-</b>
	<b>Filters</b>	
	<b>Main</b>	<b>505</b>
	<b>Reference</b>	<b>-</b>
	<b>Times</b>	
	<b>Reading 1</b>	<b>120 s</b>
	<b>Reading 2</b>	<b>-</b>
<b>Reagent 2</b>	<b>-</b>	
<b>Postdilution factor</b>	<b>2</b>	
<b>CALIBRATION</b>	<b>Type of calibration</b>	<b>Múltiplo (*)</b>
	<b>Calibrator replicates</b>	<b>3</b>
	<b>Blank replicates</b>	<b>3</b>
	<b>Calibration curve</b>	<b>-</b>
<b>OPTIONS</b>	<b>Blank absorbance limit</b>	<b>0.600</b>
	<b>Kinetic blank limit</b>	<b>-</b>
	<b>Linerity limit</b>	<b>4.5</b>

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<b>GENERAL</b>	<b>Test name</b>	<b>Microalb</b>	
	<b>Analysis mode</b>	<b>Fixed-Time</b>	
	<b>Sample type</b>	<b>Urine</b>	
	<b>Units</b>	<b>mg/L</b>	
	<b>Reaction type</b>	<b>Increasing</b>	
	<b>Turbidim. Test</b>	<b>yes</b>	
	<b>Decimals</b>	<b>0</b>	
	<b>Replicates</b>	<b>1</b>	
	<b>Name of assoc. constituent</b>	<b>-</b>	
<b>PROCEDURE</b>	<b>Volumes</b>	<b>Type of reading</b>	<b>Monoch.</b>
		<b>Sample</b>	<b>3</b>
		<b>Reagent 1</b>	<b>420</b>
		<b>Reagent 2</b>	<b>-</b>
		<b>Washing</b>	<b>1.2</b>
	<b>Filters</b>	<b>Predilution factor</b>	<b>-</b>
		<b>Main</b>	<b>535</b>
	<b>Times</b>	<b>Reference</b>	<b>-</b>
		<b>Reading 1</b>	<b>48 s</b>
		<b>Reading 2</b>	<b>168 s</b>
		<b>Reagent 2</b>	<b>-</b>
		<b>Postdilution factor</b>	<b>1.4</b>
<b>CALIBRATION</b>	<b>Type of calibration</b>	<b>specific (*)</b>	
	<b>Calibrator number</b>	<b>1</b>	
	<b>Calibrator replicates</b>	<b>3</b>	
	<b>Blank replicates</b>	<b>3</b>	
	<b>Calibration curve</b>	<b>-</b>	
<b>OPTIONS</b>	<b>Blank absorbance limit</b>	<b>1.200</b>	
	<b>Kinetic blank limit</b>	<b>-</b>	
	<b>Linerity limit</b>	<b>130</b>	

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<b>GENERAL</b>	<b>Test name</b>	<b>PCR Turbid.</b>
	<b>Analysis mode</b>	<b>Endpoint mon.</b>
	<b>Sample type</b>	<b>Serum</b>
	<b>Units</b>	<b>mg/L</b>
	<b>Reaction type</b>	<b>Increasing</b>
	<b>Turbidim. Type</b>	<b>yes</b>
	<b>Decimals</b>	<b>0</b>
	<b>Replicates</b>	<b>1</b>
	<b>Name of assoc. constituent</b>	<b>-</b>
<b>PROCEDURE</b>	<b>Type of reading</b>	<b>bichrom.</b>
	<b>Volumes</b>	
	<b>Sample</b>	<b>3</b>
	<b>Reagent 1</b>	<b>440</b>
	<b>Reagent 2</b>	<b>-</b>
	<b>Washing</b>	<b>1.2</b>
	<b>Predilution factor</b>	<b>-</b>
	<b>Filters</b>	
	<b>Main</b>	<b>535</b>
	<b>Reference</b>	<b>670</b>
	<b>Times</b>	
	<b>Reading 1</b>	<b>192 s</b>
	<b>Reading 2</b>	<b>-</b>
<b>Reagent 2</b>	<b>-</b>	
<b>Postdilution factor</b>	<b>1.4</b>	
<b>CALIBRATION</b>	<b>Type of calibration</b>	<b>specific (*)</b>
	<b>Calibrator number</b>	<b>1</b>
	<b>Calibrator replicates</b>	<b>3</b>
	<b>Blank replicates</b>	<b>3</b>
	<b>Calibration curve</b>	<b>-</b>
<b>OPTIONS</b>	<b>Blank absorbance limit</b>	<b>0.700</b>
	<b>Kinetic blank limit</b>	<b>-</b>
	<b>Linerity limit</b>	<b>150</b>

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<b>GENERAL</b>	<b>Test name</b>	<b>Proteínas totais</b>
	<b>Analysis mode</b>	<b>Endpoint mon.</b>
	<b>Sample type</b>	<b>Serum</b>
	<b>Units</b>	<b>g/dL</b>
	<b>Reaction type</b>	<b>Increasing</b>
	<b>Decimals</b>	<b>1</b>
	<b>Replicates</b>	<b>1</b>
	<b>Name of assoc. constituent</b>	<b>-</b>
<b>PROCEDURE</b>	<b>Type of reading</b>	<b>bichrom.</b>
	<b>Volumes</b>	
	<b>Sample</b>	<b>6</b>
	<b>Reagent 1</b>	<b>300</b>
	<b>Reagent 2</b>	<b>-</b>
	<b>Washing</b>	<b>1.2</b>
	<b>Predilution factor</b>	<b>-</b>
	<b>Filters</b>	
	<b>Main</b>	<b>535</b>
	<b>Reference</b>	<b>670</b>
	<b>Times</b>	
	<b>Reading 1</b>	<b>600 s</b>
	<b>Reading 2</b>	<b>-</b>
<b>Reagent 2</b>	<b>-</b>	
<b>Postdilution factor</b>	<b>2</b>	
<b>CALIBRATION</b>	<b>Type of calibration</b>	<b>Múltiplo (*)</b>
	<b>Calibrator replicates</b>	<b>3</b>
	<b>Blank replicates</b>	<b>3</b>
	<b>Calibration curve</b>	<b>-</b>
<b>OPTIONS</b>	<b>Blank absorbance limit</b>	<b>0.200</b>
	<b>Kinetic blank limit</b>	<b>-</b>
	<b>Linerity limit</b>	<b>14</b>

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<b>GENERAL</b>	<b>Test name</b>	<b>Proteinúria</b>	
	<b>Analysis mode</b>	<b>Endpoint mon.</b>	
	<b>Sample type</b>	<b>urine</b>	
	<b>Units</b>	<b>mg/dL</b>	
	<b>Reaction type</b>	<b>Increasing</b>	
	<b>Decimals</b>	<b>0</b>	
	<b>Replicates</b>	<b>1</b>	
	<b>Name of assoc. constituent</b>	<b>-</b>	
<b>PROCEDURE</b>	<b>Volumes</b>	<b>Type of reading</b>	<b>Monochrom.</b>
		<b>Sample</b>	<b>15</b>
		<b>Reagent 1</b>	<b>300</b>
		<b>Reagent 2</b>	<b>-</b>
		<b>Washing</b>	<b>1.2</b>
		<b>Predilution factor</b>	<b>-</b>
	<b>Filters</b>	<b>Main</b>	<b>600</b>
		<b>Reference</b>	
	<b>Times</b>	<b>Reading 1</b>	<b>312 s</b>
		<b>Reading 2</b>	<b>-</b>
		<b>Reagent 2</b>	<b>-</b>
		<b>Postdilution factor</b>	<b>2</b>
<b>CALIBRATION</b>	<b>Type of calibration</b>	<b>specific (*)</b>	
	<b>Calibrator replicates</b>	<b>3</b>	
	<b>Blank replicates</b>	<b>3</b>	
	<b>Calibration curve</b>	<b>-</b>	
<b>OPTIONS</b>	<b>Blank absorbance limit</b>	<b>0.300</b>	
	<b>Kinetic blank limit</b>	<b>-</b>	
	<b>Linerity limit</b>	<b>100</b>	

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<b>GENERAL</b>	<b>Test name</b>	<b>Triglicérides</b>
	<b>Analysis mode</b>	<b>Endpoint mon.</b>
	<b>Sample type</b>	<b>Serum</b>
	<b>Units</b>	<b>mg/dL</b>
	<b>Reaction type</b>	<b>Increasing</b>
	<b>Decimals</b>	<b>0</b>
	<b>Replicates</b>	<b>1</b>
	<b>Name of assoc. constituent</b>	<b>-</b>
<b>PROCEDURE</b>	<b>Type of reading</b>	<b>bichrom.</b>
	<b>Volumes</b>	
	<b>Sample</b>	<b>3</b>
	<b>Reagent 1</b>	<b>300</b>
	<b>Reagent 2</b>	<b>-</b>
	<b>Washing</b>	<b>1.2</b>
	<b>Predilution factor</b>	<b>-</b>
	<b>Filters</b>	
	<b>Main</b>	<b>505</b>
	<b>Reference</b>	<b>670</b>
	<b>Times</b>	
<b>Reading 1</b>	<b>600 s</b>	
<b>Reading 2</b>	<b>-</b>	
<b>Reagent 2</b>	<b>-</b>	
<b>Postdilution factor</b>	<b>2</b>	
<b>CALIBRATION</b>	<b>Type of calibration</b>	<b>Múltiplo (*)</b>
	<b>Calibrator replicates</b>	<b>1</b>
	<b>Blank replicates</b>	<b>1</b>
	<b>Calibration curve</b>	<b>-</b>
<b>OPTIONS</b>	<b>Blank absorbance limit</b>	<b>0.180</b>
	<b>Kinetic blank limit</b>	<b>-</b>
	<b>Linerity limit</b>	<b>1100 mg/dL</b>

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GENERAL	Test name	Uréia	
	Analysis mode	Endpoint bir.	
	Sample type	Serum	
	Units	mg/dL	
	Reaction type	Increase	
	Decimals	0	
	Replicates	1	
	Name of assoc. constituent	-	
PROCEDURE	Volumes	Type of reading	Monoch.
		Sample	3 µL
		Reagent 1	300 µL
		Reagent 2	300 µL
		Washing	1.2
		Predilution factor	-
	Filters	Main	600
		Reference	-
	Times	Reading 1	600 s
		Reading 2	-
		Reagent 2	312 s
CALIBRATION	Postdilution factor	2	
	Type of calibration	Múltiplo (*)	
	Calibrator replicates	2	
	Blank replicates	1	
	Calibration curve	-	
OPTIONS	Blank absorbance limit	0.200	
	Kinetic blank limit	-	
	Linerity limit	300 mg/dL	

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<b>GENERAL</b>	<b>Test name</b>	Uréia UV
	<b>Analysis mode</b>	Fixed-time mon.
	<b>Sample type</b>	Serum
	<b>Units</b>	mg/dL
	<b>Reaction type</b>	decreasing
	<b>Decimals</b>	0
	<b>Replicates</b>	1
	<b>Name of assoc. constituent</b>	-
<b>PROCEDURE</b>	<b>Type of reading</b>	Monoch.
	<b>Volumes</b>	
	Sample	3
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	<b>Filters</b>	
	Main	340
	Reference	-
	<b>Times</b>	
	Reading 1	48 s
Reading 2	120 s	
Reagent 2	-	
Postdilution factor	2	
<b>CALIBRATION</b>	<b>Type of calibration</b>	Múltiplo (*)
	<b>Calibrator replicates</b>	2
	<b>Blank replicates</b>	1
	<b>Calibration curve</b>	-
<b>OPTIONS</b>	<b>Blank absorbance limit</b>	0.600
	<b>Kinetic blank limit</b>	-
	<b>Linerity limit</b>	300

REVISÃO: 02/11