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Analisa

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PARAMETERS							
N°	@	Test	ACUR	Method	End-Point	Direction	Ascend
Unit	mg/dL	Decimal	1	Prim. Wavelength	510	Sec. Wavelength	700
Sample Volume	5	R1 Volume	200	R2 Volume	0	Line Limit	N/A
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	0 -- 25	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0.1	Print Number	@
Upper	@	Upper	@	Upper	25		@
Calibration							
Method	Linear	Std. 0	0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	AEO	Method	End Point	Direction	Ascend
Unit	UI/mL	Decimal	0	Prim. Wavelength	546	Sec. Wavelength	0
Sample Volume	3	R1 Volume	300	R2 Volume	0	Line Limit	N/A
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	0 -- 10	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	2	Print Number	@
Upper	@	Upper	@	Upper	800		@
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	AGPA	Method	END-POINT	Direction	Ascend
Unit	mg/dL	Decimal	0	Prim. Wavelength	340	Sec. Wavelength	0
Sample Volume	3	R1 Volume	300	R2 Volume	0	Line Limit	N/A
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	0 -- 40	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0	Print Number	@
Upper	@	Upper	@	Upper	170		
Calibration							
Method	Curve	Std. 0	0				
		Std. 1	#				
		Std. 2	#				
		Std. 3	#				
		Std. 4	#				
		Std. 5	#				

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PARAMETERS							
N°	@	Test	ALB	Method	End-Point	Direction	Ascend
Unit	g/dL	Decimal	1	Prim. Wavelength	630	Sec. Wavelength	0
Sample Volume	3	R1 Volume	300	R2 Volume	0	Line Limit	N/A
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	0 -- 10	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0.1	Print Number	@
Upper	@	Upper	@	Upper	7.0		@
Calibration							
Method	Linear	Std. 0	0				
		Std. 1	#				
		Std. 2	#				
		Std. 3	#				
		Std. 4	#				
		Std. 5	#				

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PARAMETERS							
N°	@	Test	ALT	Method	Kinetic	Direction	Descend
Unit	U/L	Decimal	0	Prim. Wavelength	340	Sec. Wavelength	0
Sample Volume	10	R1 Volume	200	R2 Volume	0	Line Limit	20
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	5 -- 20	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	1.6	Print Number	@
Upper	@	Upper	@	Upper	500		@
Calibration							
Method	Linear	Std. 0	0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	AMIL	Method	Kinetic	Direction	Ascend
Unit	U/L	Decimal	0	Prim. Wavelength	405	Sec. Wavelength	0
Sample Volume	4	R1 Volume	200	R2 Volume	0	Line Limit	20
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	5 -- 20	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	2	Print Number	@
Upper	@	Upper	@	Upper	1300		@
Calibration							
Method	Linear	Std. 0	0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	AST	Method	Kinetic	Direction	Descend
Unit	U/L	Decimal	0	Prim. Wavelength	340	Sec. Wavelength	0
Sample Volume	10	R1 Volume	200	R2 Volume	0	Line Limit	20
Incubation	--	Check Antigen	No	Response		Sample Volume	@
Reaction	5 -- 20	Substrate	N/A	Lower	--	Dilution	@
				Upper	--		
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	1.6	Print Number	@
Upper	@	Upper	@	Upper	500		@
Calibration							
Method	Linear	Std. 0	0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	BILT	Method	End Point	Direction	Ascend
Unit	mg/dL	Decimal	1	Prim. Wavelength	546	Sec. Wavelength	700
Sample Volume	20	R1 Volume	200	R2 Volume	0	Line Limit	N/A
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	0 -- 10	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0.02	Print Number	@
Upper	@	Upper	@	Upper	15		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	Ca Ars	Method	End Point	Direction	Ascend
Unit	mg/dL	Decimal	1	Prim. Wavelength	630	Sec. Wavelength	0
Sample Volume	4	R1 Volume	250	R2 Volume	0	Line Limit	N/A
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	0 -- 25	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0.2	Print Number	@
Upper	@	Upper	@	Upper	18		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	Ca MTB	Method	End Point	Direction	Ascend
Unit	mg/dL	Decimal	1	Prim. Wavelength	630	Sec. Wavelength	0
Sample Volume	2	R1 Volume	200	R2 Volume	0	Line Limit	N/A
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	0 -- 10	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity			
Lower	@	Lower	@	Lower	0.6	Full Name	@
Upper	@	Upper	@	Upper	15	Print Number	@
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	COL	Method	End Point	Direction	Ascend
Unit	mg/dL	Decimal	0	Prim. Wavelength	510	Sec. Wavelength	0
Sample Volume	3	R1 Volume	300	R2 Volume	0	Line Limit	N/A
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	0 -- 25	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity			
Lower	@	Lower	@	Lower	1	Full Name	@
Upper	@	Upper	@	Upper	1000	Print Number	@
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	HDL	Method	End Point	Direction	Ascend
Unit	mg/dL	Decimal	0	Prim. Wavelength	510	Sec. Wavelength	0
Sample Volume	20	R1 Volume	200	R2 Volume	0	Line Limit	N/A
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	0 -- 50	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	2	Print Number	@
Upper	@	Upper	@	Upper	150		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	CHE	Method	Kinetic	Direction	Descend
Unit	U/L	Decimal	0	Prim. Wavelength	405	Sec. Wavelength	0
Sample Volume	5	R1 Volume	300	R2 Volume	0	Line Limit	20
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	5 -- 20	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	123	Print Number	@
Upper	@	Upper	@	Upper	25000		
Calibration							
Method	Factor						
Factor	65804						

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PARAMETERS							
N°	@	Test	CK	Method	Kinetic	Direction	Ascend
Unit	U/L	Decimal	0	Prim. Wavelength	340	Sec. Wavelength	0
Sample Volume	10	R1 Volume	200	R2 Volume	0	Line Limit	20
Incubation	--	Check Antigen	No	Response		Sample Volume	@
Reaction	15 -- 30	Substrate	N/A	Lower	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	10	Print Number	@
Upper	@	Upper	@	Upper	900		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	CK MB	Method	Fixed Time	Direction	Ascend
Unit	U/L	Decimal	0	Prim. Wavelength	340	Sec. Wavelength	0
Sample Volume	8	R1 Volume	200	R2 Volume	0	Line Limit	20
Incubation	--	Check Antigen	No	Response		Sample Volume	@
Reaction	50 -- 75	Substrate	N/A	Lower	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	3	Print Number	@
Upper	@	Upper	@	Upper	330		
Calibration							
Method	Factor						
Factor	1650						

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PARAMETERS							
N°	@	Test	CREA	Method	Fixed Time	Direction	Ascend
Unit	mg/dL	Decimal	1	Prim. Wavelength	510	Sec. Wavelength	0
Sample Volume	20	R1 Volume	200	R2 Volume	0	Line Limit	20
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	3 -- 8	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0.1	Print Number	@
Upper	@	Upper	@	Upper	20		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	FERR	Method	Fixed Time	Direction	Ascend
Unit	µg/L	Decimal	0	Prim. Wavelength	546	Sec. Wavelength	0
Sample Volume	6	R1 Volume	200	R2 Volume	0	Line Limit	20
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	3 -- 28	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0	Print Number	@
Upper	@	Upper	@	Upper	500		
Calibration							
Method	Curve	Std. 0	0				
		Std. 1	#				
		Std. 2	#				
		Std. 3	#				
		Std. 4	#				
		Std. 5	#				

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PARAMETERS							
N°	@	Test	FERC	Method	End Point	Direction	Ascend
Unit	µg/dL	Decimal	0	Prim. Wavelength	630	Sec. Wavelength	0
Sample Volume	15	R1 Volume	300	R2 Volume	0	Line Limit	N/A
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	0 -- 50	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	10	Print Number	@
Upper	@	Upper	@	Upper	500		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	FERF	Method	End Point	Direction	Ascend
Unit	µg/dL	Decimal	1	Prim. Wavelength	578	Sec. Wavelength	0
Sample Volume	40	R1 Volume	160	R2 Volume	40	Line Limit	N/A
				Response			
Incubation	5	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	-1 -- 25	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	10	Print Number	@
Upper	@	Upper	@	Upper	1000		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	FALC	Method	Kinetic	Direction	Ascend
Unit	U/L	Decimal	0	Prim. Wavelength	405	Sec. Wavelength	0
Sample Volume	4	R1 Volume	200	R2 Volume	0	Line Limit	20
Incubation	--	Check Antigen	No	Response		Sample Volume	@
Reaction	5 -- 20	Substrate	N/A	Lower	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	1	Print Number	@
Upper	@	Upper	@	Upper	1200		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	FALC	Method	Kinetic	Direction	Ascend
Unit	U/L	Decimal	0	Prim. Wavelength	405	Sec. Wavelength	0
Sample Volume	4	R1 Volume	200	R2 Volume	0	Line Limit	20
Incubation	--	Check Antigen	No	Response		Sample Volume	@
Reaction	5 -- 20	Substrate	N/A	Lower	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	1	Print Number	@
Upper	@	Upper	@	Upper	690		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	FOSF	Method	End Point	Direction	Ascend
Unit	mg/dL	Decimal	1	Prim. Wavelength	340	Sec. Wavelength	0
Sample Volume	3	R1 Volume	210	R2 Volume	90	Line Limit	N/A
				Response			
Incubation	10	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	-1 -- 25	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0.1	Print Number	@
Upper	@	Upper	@	Upper	20		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	FR	Method	End-Point	Direction	Ascend
Unit	UI/mL	Decimal	0	Prim. Wavelength	630	Sec. Wavelength	0
Sample Volume	3	R1 Volume	240	R2 Volume	60	Line Limit	N/A
				Response			
Incubation	5	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	-1 -- 10	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	2	Print Number	@
Upper	@	Upper	@	Upper	160		
Calibration							
Method	Curve	Std. 0	0				
		Std. 1	#				
		Std. 2	#				
		Std. 3	#				
		Std. 4	#				
		Std. 5	#				

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PARAMETERS							
N°	@	Test	FRUT	Method	Fixed Time	Direction	Ascend
Unit	mmol/L	Decimal	1	Prim. Wavelength	546	Sec. Wavelength	0
Sample Volume	10	R1 Volume	200	R2 Volume	0	Line Limit	20
Incubation	--	Check Antigen	No	Response	Lower --	Sample Volume	@
Reaction	50 -- 75	Substrate	N/A		Upper --	Dilution	@
R1 Blank		Mix R Blank		Linearity	Lower 0.1	Full Name	@
	Lower @	Lower @			Upper 7.0	Print Number	@
	Upper @	Upper @					
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	GGT	Method	Kinetic	Direction	Ascend
Unit	U/L	Decimal	0	Prim. Wavelength	405	Sec. Wavelength	0
Sample Volume	20	R1 Volume	200	R2 Volume	0	Line Limit	20
Incubation	--	Check Antigen	No	Response		Sample Volume	@
Reaction	5 -- 20	Substrate	N/A	Lower	--	Dilution	@
R1 Blank		Mix R Blank		Upper	--	Full Name	@
Lower	@	Lower	@	Linearity		Print Number	@
Upper	@	Upper	@	Lower	2		
Calibration				Upper	300		
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	GLIC	Method	End Point	Direction	Ascend
Unit	mg/dL	Decimal	0	Prim. Wavelength	510	Sec. Wavelength	0
Sample Volume	3	R1 Volume	300	R2 Volume	0	Line Limit	N/A
Incubation	--	Check Antigen	No	Response		Sample Volume	@
Reaction	0 -- 25	Substrate	N/A	Lower	--	Dilution	@
R1 Blank		Mix R Blank		Upper	--	Full Name	@
Lower	@	Lower	@	Linearity		Print Number	@
Upper	@	Upper	@	Lower	1		
Calibration				Upper	500		
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	HDL	Method	End Point	Direction	Ascend
Unit	mg/dL	Decimal	0	Prim. Wavelength	578	Sec. Wavelength	700
Sample Volume	2	R1 Volume	200	R2 Volume	67	Line Limit	N/A
Response							
Incubation	25	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	-1 -- 25	Substrate	N/A	Upper	--	Dilution	@
Linearity							
R1 Blank		Mix R Blank		Lower	1	Full Name	@
Lower	@	Lower	@	Upper	200	Print Number	@
Upper	@	Upper	@				
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	LDL	Method	End Point	Direction	Ascend
Unit	mg/dL	Decimal	0	Prim. Wavelength	546	Sec. Wavelength	700
Sample Volume	2	R1 Volume	200	R2 Volume	67	Line Limit	N/A
Response							
Incubation	25	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	-1 -- 25	Substrate	N/A	Upper	--	Dilution	@
Linearity							
R1 Blank		Mix R Blank		Lower	1	Full Name	@
Lower	@	Lower	@	Upper	990	Print Number	@
Upper	@	Upper	@				
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	LDH	Method	Kinetic	Direction	Descend
Unit	U/L	Decimal	0	Prim. Wavelength	340	Sec. Wavelength	0
Sample Volume	4	R1 Volume	200	R2 Volume	0	Line Limit	20
Response							
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	5 -- 20	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	2	Print Number	@
Upper	@	Upper	@	Upper	1250		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	Mg	Method	End Point	Direction	Ascend
Unit	mg/dL	Decimal	1	Prim. Wavelength	510	Sec. Wavelength	0
Sample Volume	3	R1 Volume	300	R2 Volume	0	Line Limit	N/A
Response							
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	0 -- 10	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0.1	Print Number	@
Upper	@	Upper	@	Upper	4.0		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test		Method	Fixed-Time	Direction	Ascend
Unit	mg/L	Decimal	1	Prim. Wavelength	546	Sec. Wavelength	0
Sample Volume	2	R1 Volume	286	R2 Volume	0	Line Limit	20
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	3 -- 14	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	1	Print Number	@
Upper	@	Upper	@	Upper	130		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	PCR	Method	End Point	Direction	Ascend
Unit	mg/L	Decimal	1	Prim. Wavelength	546	Sec. Wavelength	0
Sample Volume	2	R1 Volume	286	R2 Volume	0	Line Limit	N/A
Incubation	--	Check Antigen	No	Response		Sample Volume	@
Reaction	0 -- 15	Substrate	N/A	Lower	--	Dilution	@
R1 Blank		Mix R Blank		Upper	--		
Lower	@	Lower	@	Linearity		Full Name	@
Upper	@	Upper	@	Lower	1	Print Number	@
				Upper	150		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	PCR-AS	Method	Fixed-time	Direction	Ascend
Unit	mg/L	Decimal	2	Prim. Wavelength	546	Sec. Wavelength	0
Sample Volume	3	R1 Volume	225	R2 Volume	0	Line Limit	20
Incubation	--	Check Antigen	No	Response		Sample Volume	@
Reaction	3 -- 28	Substrate	N/A	Lower	--	Dilution	@
R1 Blank		Mix R Blank		Upper	--		
Lower	@	Lower	@	Linearity		Full Name	@
Upper	@	Upper	@	Lower	0.06	Print Number	@
				Upper	15		
Calibration							
Method	Curve	Std. 0	0.0				
		Std. 1	#				
		Std. 2	#				
		Std. 3	#				
		Std. 4	#				
		Std. 5	#				

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PARAMETERS							
N°	@	Test	PROT	Method	End Point	Direction	Ascend
Unit	g/dL	Decimal	1	Prim. Wavelength	546	Sec. Wavelength	0
Sample Volume	4	R1 Volume	200	R2 Volume	0	Line Limit	N/A
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	0 -- 25	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0.02	Print Number	@
Upper	@	Upper	@	Upper	15		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	PTur	Method	End Point	Direction	Ascend
Unit	mg/dL	Decimal	0	Prim. Wavelength	630	Sec. Wavelength	0
Sample Volume	10	R1 Volume	200	R2 Volume	0	Line Limit	N/A
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	0 -- 50	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	@
Lower	@	Lower	@	Lower	0.02	Print Number	@
Upper	@	Upper	@	Upper	160		
Calibration							
Method	Linear	Std. 0	0.0				
		Std. 1	#				
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	TRI	Method	End Point	Direction	Ascend
Unit	mg/dL	Decimal	0	Prim. Wavelength	510	Sec. Wavelength	700
Sample Volume	2	R1 Volume	200	R2 Volume	0	Line Limit	N/A
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	0 -- 25	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	
Lower	@	Lower	@	Lower	1	Print Number	@
Upper	@	Upper	@	Upper	600		
Calibration Method		Std. 0 0.0					
		Std. 1 #					
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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PARAMETERS							
N°	@	Test	URE	Method	Fixed Time	Direction	Descend
Unit	mg/dL	Decimal	1	Prim. Wavelength	340	Sec. Wavelength	0
Sample Volume	2	R1 Volume	300	R2 Volume	0	Line Limit	20
				Response			
Incubation	--	Check Antigen	No	Lower	--	Sample Volume	@
Reaction	3 -- 8	Substrate	N/A	Upper	--	Dilution	@
R1 Blank		Mix R Blank		Linearity		Full Name	
Lower	@	Lower	@	Lower	1	Print Number	@
Upper	@	Upper	@	Upper	300		
Calibration Method		Std. 0 0.0					
		Std. 1 #					
		Std. 2					
		Std. 3					
		Std. 4					
		Std. 5					

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