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## Instrument Settings

Test Methodology	
Kind of Process	<b>End Point</b>
Type	<b>Linear</b>
Filters	<b>510 / 630</b>
Reaction Direction	<b>Increase</b>
Reagent #1 (µL)	<b>300</b>
Reagent #2 (µL)	-
Sample Starter	<b>Inactive</b>
Delay Time (sec.)	<b>0</b>
Incubation Time (sec.)	<b>300</b>
Reading Time (sec.)	<b>20</b>

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## Instrument Settings

Test Methodology	<b>Turbidimetry</b>
Kind of Process	<b>Fixed time</b>
Type	<b>Linear</b>
Filters	<b>546</b>
Reaction Direction	<b>Increase</b>
Reagent #1 (µL)	<b>300</b>
Reagent #2 (µL)	<b>-</b>
Sample Starter	<b>Inactive</b>
Delay Time (sec.)	<b>10</b>
Incubation Time (sec.)	<b>30 - 0</b>
Reading Time (sec.)	<b>120</b>

Units – Serum	<b>UI/mL</b>
Units – Urine	
Number of needle washes	<b>2/2</b>
Number of cuvette washes	<b>2</b>
Dynamic Blank	<b>Inactive</b>
Reagent Blank	<b>Every Run</b>
Reagent Limit (mABS)	<b>1500</b>
Curve Acceptance (%)	<b>100</b>
Instrument Factor	<b>0.000</b>
Shift	<b>0.000</b>

<b>SERUM</b>	
Name	<b>AEO</b>
Sample (µL)	<b>3</b>
Pre-Dilution	<b>1.00</b>
<u>Dilution</u>	
Factor	<b>2.00</b>
Test Limit	<b>800</b>
Initial ABS	<b>0</b>
Final ABS	<b>1900</b>
Max ABS Delta	<b>500</b>
<u>Normal Values</u>	
Man	<b>0 – 200</b>
Woman	<b>0 – 200</b>
Child	<b>0 – 150</b>
Calibrator	
Std. 1	<b>0</b>
Std. 2	<b>#</b>

<b>URINE</b>	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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## Instrument Settings

Test Methodology	<b>DCA</b>
Kind of Process	<b>End Point 2 points</b>
Type	<b>Linear</b>
Filters	<b>546/630</b>
Reaction Direction	<b>Increase</b>
Reagent #1 (µL)	<b>300</b>
Reagent #2 (µL)	<b>75</b>
Sample Starter	<b>Inactive</b>
Delay Time (sec.)	<b>0</b>
Incubation Time (sec.)	<b>20 - 20</b>
Reading Time (sec.)	<b>180</b>

Units – Serum	<b>mg/dL</b>
Units – Urine	<b>---</b>
Number of needle washes	<b>1/1</b>
Number of cuvette washes	<b>1</b>
Dynamic Blank	<b>Inactive</b>
Reagent Blank	<b>Every Run</b>
Reagent Limit (mABS)	<b>300</b>
Curve Acceptance (%)	<b>20</b>
Instrument Factor	<b>0.000</b>
Shift	<b>0.000</b>

<b>SERUM</b>	
Name	<b>BIL-D</b>
Sample (µL)	<b>30</b>
Pré-Dilution	<b>1.00</b>
<u>Dilution</u>	
Factor	<b>10.00</b>
Test Limit	<b>12</b>
Initial ABS	<b>@</b>
Final ABS	<b>@</b>
Max ABS Delta	<b>99999</b>
<u>Normal Values</u>	
Man	<b>0 – 0.2</b>
Woman	<b>0 – 0.2</b>
Child	<b>-</b>
Calibrator	
Std. 1	<b>0</b>
Std. 2	<b>#</b>

<b>URINE</b>	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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## Instrument Settings

Test Methodology	
Kind of Process	<b>End Point 2 points</b>
Type	<b>Linear</b>
Filters	<b>546/630</b>
Reaction Direction	<b>Increase</b>
Reagent #1 (µL)	<b>300</b>
Reagent #2 (µL)	<b>75</b>
Sample Starter	<b>Inactive</b>
Delay Time (sec.)	<b>0</b>
Incubation Time (sec.)	<b>30 - 240</b>
Reading Time (sec.)	<b>20</b>

Units – Serum	<b>mg/dL</b>
Units – Urine	
Number of needle washes	<b>1/1</b>
Number of cuvette washes	<b>1</b>
Dynamic Blank	<b>Inactive</b>
Reagent Blank	<b>Every Run</b>
Reagent Limit (mABS)	<b>300</b>
Curve Acceptance (%)	<b>20</b>
Instrument Factor	<b>0.000</b>
Shift	<b>0.000</b>

SERUM	
Name	<b>BIL-T</b>
Sample (µL)	<b>19</b>
Pré-Dilution	<b>1.00</b>
<u>Dilution</u>	
Factor	<b>10.00</b>
Test Limit	<b>30</b>
Initial ABS	<b>@</b>
Final ABS	<b>@</b>
Max ABS Delta	<b>99999</b>
<u>Normal Values</u>	
Man	<b>0 – 1.2</b>
Woman	<b>0 – 1.2</b>
Child	<b>-</b>
Calibrator	
Std. 1	<b>0</b>
Std. 2	<b>#</b>

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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## Instrument Settings

Test Methodology	
Kind of Process	<b>End Point</b>
Type	<b>Linear</b>
Filters	<b>510/630</b>
Reaction Direction	<b>Increase</b>
Reagent #1 (µL)	<b>300</b>
Reagent #2 (µL)	<b>-</b>
Sample Starter	<b>Inactive</b>
Delay Time (sec.)	<b>0</b>
Incubation Time (sec.)	<b>300</b>
Reading Time (sec.)	<b>20</b>

Units – Serum	<b>mg/dL</b>
Units – Urine	
Number of needle washes	<b>1/1</b>
Number of cuvette washes	<b>1</b>
Dynamic Blank	<b>Inactive</b>
Reagent Blank	<b>Every Run</b>
Reagent Limit (mABS)	<b>250</b>
Curve Acceptance (%)	<b>20</b>
Instrument Factor	<b>0.000</b>
Shift	<b>0.000</b>

SERUM	
Name	<b>CHOL</b>
Sample (µL)	<b>3</b>
Pré-Dilution	<b>1.00</b>
<u>Dilution</u>	
Factor	<b>2.00</b>
Test Limit	<b>500</b>
Initial ABS	<b>@</b>
Final ABS	<b>@</b>
Max ABS Delta	<b>99999</b>
<u>Normal Values</u>	
Man	<b>&lt; 200</b>
Woman	<b>&lt; 200</b>
Child	<b>-----</b>
Calibrator	
Std. 1	<b>0</b>
Std. 2	<b>#</b>

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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### Instrument Settings

Test Methodology	
Kind of Process	<b>Fixed Time</b>
Type	<b>Linear</b>
Filters	<b>510</b>
Reaction Direction	<b>Increase</b>
Reagent #1 (µL)	<b>60</b>
Reagent #2 (µL)	<b>240</b>
Sample Starter	<b>Inactive</b>
Delay Time (sec.)	<b>0</b>
Incubation Time (sec.)	<b>60 - 0</b>
Reading Time (sec.)	<b>120</b>

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## Instrument Settings

Test Methodology	
Kind of Process	<b>End Point 2 points</b>
Type	<b>Linear</b>
Filters	<b>578</b>
Reaction Direction	<b>Increase</b>
Reagent #1 (µL)	<b>240</b>
Reagent #2 (µL)	<b>60</b>
Sample Starter	<b>Inactive</b>
Delay Time (sec.)	<b>0</b>
Incubation Time (sec.)	<b>60 - 0</b>
Reading Time (sec.)	<b>300</b>

Units – Serum	<b>µg/dL</b>
Units – Urine	
Number of needle washes	<b>2/2</b>
Number of cuvette washes	<b>2</b>
Dynamic Blank	<b>Inactive</b>
Reagent Blank	<b>Every Run</b>
Reagent Limit (mABS)	<b>1500</b>
Curve Acceptance (%)	<b>100</b>
Instrument Factor	<b>0.000</b>
Shift	<b>0.000</b>

SERUM	
Name	<b>Fe</b>
Sample (µL)	<b>30</b>
Pré-Dilution	<b>1.00</b>
<u>Dilution</u>	
Factor	<b>2.00</b>
Test Limit	<b>1000</b>
Initial ABS	<b>1900</b>
Final ABS	<b>1900</b>
Max ABS Delta	<b>500</b>
<u>Normal Values</u>	
Man	<b>65 - 170</b>
Woman	<b>50 - 170</b>
Child	<b>-----</b>
Calibrator	
Std. 1	<b>0</b>
Std. 2	<b>#</b>

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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## Instrument Settings

Test Methodology	
Kind of Process	<b>Kinetic</b>
Type	<b>Linear</b>
Filters	<b>405</b>
Reaction Direction	<b>Increase</b>
Reagent #1 (µL)	<b>300</b>
Reagent #2 (µL)	<b>----</b>
Sample Starter	<b>Inactive</b>
Delay Time (sec.)	<b>0</b>
Incubation Time (sec.)	<b>60</b>
Reading Time (sec.)	<b>180</b>

Units – Serum	<b>U/L</b>
Units – Urine	
Number of needle washes	<b>2/2</b>
Number of cuvette washes	<b>2</b>
Dynamic Blank	<b>Inactive</b>
Reagent Blank	<b>Every Run</b>
Reagent Limit (mABS)	<b>1500</b>
Curve Acceptance (%)	<b>100</b>
Instrument Factor	<b>0.000</b>
Shift	<b>0.000</b>

<b>SERUM</b>	
Name	<b>FALC</b>
Sample (µL)	<b>6</b>
Pré-Dilution	<b>1.00</b>
<u>Dilution</u>	
Factor	<b>2.00</b>
Test Limit	<b>1500</b>
Initial ABS	<b>@</b>
Final ABS	<b>@</b>
Max ABS Delta	<b>500</b>
<u>Normal Values</u>	
Man	<b>27 – 100</b>
Woman	<b>27 – 100</b>
Child	<b>-</b>
Calibrator	
Std. 1	<b>0</b>
Std. 2	<b>#</b>

<b>URINE</b>	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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Test Methodology	<b>Molibd UV</b>
Method	<b>End Point</b>
Kind of Process	<b>Linear</b>
Filters	<b>340/700</b>
Reaction Direction	<b>Increase</b>
Reagent #1 (µL)	<b>300</b>
Reagent #2 (µL)	<b>0</b>
Sample Starter	<b>(Inactive)</b>
Delay Time (sec.)	<b>0</b>
Incubation Time (sec.)	<b>300</b>
Reading Time (sec.)	<b>10</b>

Units – Serum	<b>mg/dL</b>
Units – Urine	<b>mg/dL</b>
Number of needle washes	<b>2</b>
Number of cuvette washes	<b>2</b>
Dynamic Blank	<b>Inactive</b>
Reagent Blank	<b>Every Run</b>
Reagent Limit (mABS)	<b>1500</b>
Curve Acceptance (%)	<b>100</b>

<b>SERUM</b>	
Name	<b>Fósforo</b>
Sample (µL)	<b>3</b>
Pré-Dilution	<b>1.00</b>
<u>Dilution</u>	
Factor	<b>2.00</b>
Test Limit (conc)	<b>20</b>
Max ABS Delta (mABS)	<b>2000</b>
Check Prozone	<b>inactive</b>
Instrumental Factor	<b>0.000</b>
Shift	<b>0.000</b>
Re-Run Hyperactive	<b>active</b>
Re-Run Pathological	<b>Inactive</b>
<u>Normal range</u>	
Male	<b>2.50 – 4.80</b>
Female	<b>2.50 – 4.80</b>
Child	<b>4.50 – 6.70</b>

<b>URINE</b>	
Name	<b>Fósforo</b>
Sample (µL)	<b>3</b>
Pré-Dilution	<b>20.00</b>
<u>Dilution</u>	
Factor	<b>1.00</b>
Test Limit (conc)	<b>20</b>
Check Prozone	<b>inactive</b>
Max ABS Delta (mABS)	<b>2000</b>
Instrumental Factor	<b>0.000</b>
Shift	<b>0.000</b>
Re-Run Hyperactive	<b>active</b>
Re-Run Pathological	<b>Inactive</b>
<u>Normal range</u>	
Male	<b>0.00 / 0.00</b>
Female	<b>0.00 / 0.00</b>
Child	<b>0.00 / 0.00</b>

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## Instrument Settings

Test Methodology	
Kind of Process	<b>Fixed Time</b>
Type	<b>Linear</b>
Filters	<b>546</b>
Reaction Direction	<b>Increase</b>
Reagent #1 (µL)	<b>300</b>
Reagent #2 (µL)	<b>----</b>
Sample Starter	<b>Inactive</b>
Delay Time (sec.)	<b>0</b>
Incubation Time (sec.)	<b>600</b>
Reading Time (sec.)	<b>300</b>

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## Instrument Settings

Test Methodology	
Kind of Process	<b>End Point</b>
Type	<b>Linear</b>
Filters	<b>510</b>
Reaction Direction	<b>Increase</b>
Reagent #1 (µL)	<b>300</b>
Reagent #2 (µL)	<b>-</b>
Sample Starter	<b>Inactive</b>
Delay Time (sec.)	<b>0</b>
Incubation Time (sec.)	<b>300</b>
Reading Time (sec.)	<b>10</b>

Units – Serum	<b>mg/dL</b>
Units – Urine	
Number of needle washes	<b>1</b>
Number of cuvette washes	<b>1</b>
Dynamic Blank	<b>Inactive</b>
Reagent Blank	<b>Every Run</b>
Reagent Limit (mABS)	<b>250</b>
Curve Acceptance (%)	<b>20</b>
Instrument Factor	<b>0.000</b>
Shift	<b>0.000</b>

SERUM	
Name	<b>GLU</b>
Sample (µL)	<b>3</b>
Pré-Dilution	<b>1.00</b>
<u>Dilution</u>	
Factor	<b>2.00</b>
Test Limit	<b>500</b>
Initial ABS	<b>@</b>
Final ABS	<b>@</b>
Max ABS Delta	<b>500</b>
<u>Normal Values</u>	
Man	<b>70 - 105</b>
Woman	<b>70 - 105</b>
Child	<b>70 - 90</b>
Calibrator	
Std. 1	<b>0</b>
Std. 2	<b>#</b>

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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Test Methodology	<b>Turbidimetry</b>
Kind of Process	<b>End Point</b>
Type	<b>Linear</b>
Filters	<b>546/700</b>
Reaction Direction	<b>Increase</b>
Reagent #1 (µL)	<b>300</b>
Reagent #2 (µL)	<b>-</b>
Sample Starter	<b>Inactive</b>
Delay Time (sec.)	<b>0</b>
Incubation Time (sec.)	<b>300</b>
Reading Time (sec.)	<b>10</b>

Units – Serum	<b>mmol/L</b>
Units – Urine	
Number of needle washes	<b>2/2</b>
Number of cuvette washes	<b>2</b>
Dynamic Blank	<b>Inactive</b>
Reagent Blank	<b>Every Run</b>
Reagent Limit (mABS)	<b>1500</b>
Curve Acceptance (%)	<b>100</b>
Instrument Factor	<b>0.000</b>
Shift	<b>0.000</b>

<b>SERUM</b>	
Name	<b>Hb</b>
Sample (µL)	<b>50</b>
Pré-Dilution	<b>1.00</b>
<u>Dilution</u>	
Factor	<b>2.00</b>
Test Limit	<b>24</b>
Initial ABS	<b>0</b>
Final ABS	<b>1900</b>
Max ABS Delta	<b>500</b>
<u>Normal Values</u>	
Man	
Woman	
Child	
Calibrator	
Std. 1	<b>0.000</b>
Std. 2	<b>#</b>

<b>URINE</b>	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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Test Methodology	<b>Turbidimetry</b>
Kind of Process	<b>End Point 2 points</b>
Type	<b>Cubic Spline</b>
Filters	<b>340</b>
Reaction Direction	<b>Increase</b>
Reagent #1 (µL)	<b>250</b>
Reagent #2 (µL)	<b>50</b>
Sample Starter	<b>Inactive</b>
Delay Time (sec.)	<b>10</b>
Incubation Time (sec.)	<b>300 - 0</b>
Reading Time (sec.)	<b>660</b>

Units – Serum	<b>mmol/L</b>
Units – Urine	
Number of needle washes	<b>2/2</b>
Number of cuvette washes	<b>2</b>
Dynamic Blank	<b>Inactive</b>
Reagent Blank	<b>Every Run</b>
Reagent Limit (mABS)	<b>1600</b>
Curve Acceptance (%)	<b>100</b>
Instrument Factor	<b>0.000</b>
Shift	<b>0.000</b>

SERUM	
Name	<b>HbA1c</b>
Sample (µL)	<b>10</b>
Pré-Dilution	<b>1.00</b>
<u>Dilution</u>	
Factor	<b>2.00</b>
Test Limit	<b>2.50</b>
Initial ABS	<b>1900</b>
Final ABS	<b>1900</b>
Max ABS Delta	<b>500</b>
<u>Normal Values</u>	
Man	
Woman	
Child	
Calibrator	
Std. 1	<b>0.000</b>
Std. 2	<b>#</b>
Std. 3	<b>#</b>
Std. 4	<b>#</b>
Std.5	<b>#</b>
Std.6	<b>#</b>

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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## Instrument Settings

Test Methodology	
Kind of Process	<b>End Point 2 points</b>
Type	<b>Linear</b>
Filters	<b>578/700</b>
Reaction Direction	<b>Increase</b>
Reagent #1 (µL)	<b>300</b>
Reagent #2 (µL)	<b>100</b>
Sample Starter	<b>Inactive</b>
Delay Time (sec.)	<b>0</b>
Incubation Time (sec.)	<b>300 - 0</b>
Reading Time (sec.)	<b>300</b>

Units – Serum	<b>mg/dL</b>
Units – Urine	
Number of needle washes	<b>2/2</b>
Number of cuvette washes	<b>2</b>
Dynamic Blank	<b>Inactive</b>
Reagent Blank	<b>Every Run</b>
Reagent Limit (mABS)	<b>1500</b>
Curve Acceptance (%)	<b>100</b>
Instrument Factor	<b>0.000</b>
Shift	<b>0.000</b>

SERUM	
Name	<b>HDL-D</b>
Sample (µL)	<b>3</b>
Pré-Dilution	<b>1.00</b>
<u>Dilution</u>	
Factor	<b>2.00</b>
Test Limit	<b>200</b>
Initial ABS	<b>@</b>
Final ABS	<b>@</b>
Max ABS Delta	<b>500</b>
<u>Normal Values</u>	
Man	<b>&gt; 55</b>
Woman	<b>&gt; 65</b>
Child	<b>-----</b>
Calibrator	
Std. 1	<b>0</b>
Std. 2	<b>#</b>

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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### Instrument Settings

Test Methodology	
Kind of Process	<b>End Point 2 points</b>
Type	<b>Linear</b>
Filters	<b>546</b>
Reaction Direction	<b>Increase</b>
Reagent #1 (µL)	<b>300</b>
Reagent #2 (µL)	<b>100</b>
Sample Starter	<b>Inactive</b>
Delay Time (sec.)	<b>0</b>
Incubation Time (sec.)	<b>300 - 0</b>
Reading Time (sec.)	<b>300</b>

Units – Serum	<b>mg/dL</b>
Units – Urine	
Number of needle washes	<b>2/2</b>
Number of cuvette washes	<b>2</b>
Dynamic Blank	<b>Inactive</b>
Reagent Blank	<b>Every Run</b>
Reagent Limit (mABS)	<b>1500</b>
Curve Acceptance (%)	<b>100</b>
Instrument Factor	<b>0.000</b>
Shift	<b>0.000</b>

SERUM	
Name	<b>LDL-D</b>
Sample (µL)	<b>3</b>
Pré-Dilution	<b>1.00</b>
<u>Dilution</u>	
Factor	<b>2.00</b>
Test Limit	<b>990</b>
Initial ABS	<b>@</b>
Final ABS	<b>@</b>
Max ABS Delta	<b>500</b>
<u>Normal Values</u>	
Man	<b>&lt; 130</b>
Woman	<b>&lt; 130</b>
Child	<b>-----</b>
Calibrator	
Std. 1	<b>0</b>
Std. 2	<b>#</b>

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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## Instrument Settings

Test Methodology	
Kind of Process	<b>End Point</b>
Type	<b>Linear</b>
Filters	<b>510</b>
Reaction Direction	<b>Increase</b>
Reagent #1 (µL)	<b>300</b>
Reagent #2 (µL)	<b>-</b>
Sample Starter	<b>Inactive</b>
Delay Time (sec.)	<b>10</b>
Incubation Time (sec.)	<b>120</b>
Reading Time (sec.)	<b>20</b>

Units – Serum	<b>mg/dL</b>
Units – Urine	
Number of needle washes	<b>2</b>
Number of cuvette washes	<b>2</b>
Dynamic Blank	<b>Inactive</b>
Reagent Blank	<b>Every Run</b>
Reagent Limit (mABS)	<b>800</b>
Curve Acceptance (%)	<b>100</b>
Instrument Factor	<b>0.000</b>
Shift	<b>0.000</b>

SERUM	
Name	<b>MG</b>
Sample (µL)	<b>3</b>
Pré-Dilution	<b>1.00</b>
<u>Dilution</u>	
Factor	<b>2.00</b>
Test Limit	<b>4.5</b>
Initial ABS	<b>@</b>
Final ABS	<b>@</b>
Max ABS Delta	<b>99999</b>
<u>Normal Values</u>	
Man	<b>1.9 – 2.5</b>
Woman	<b>1.9 – 2.5</b>
Child	<b>@</b>
Calibrator	
Std. 1	<b>0</b>
Std. 2	<b>#</b>

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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## Instrument Settings

Test Methodology	
Kind of Process	<b>End Point</b>
Type	<b>Linear</b>
Filters	<b>546</b>
Reaction Direction	<b>Increase</b>
Reagent #1 (µL)	<b>300</b>
Reagent #2 (µL)	<b>-</b>
Sample Starter	<b>Inactive</b>
Delay Time (sec.)	<b>0</b>
Incubation Time (sec.)	<b>600</b>
Reading Time (sec.)	<b>10</b>

Units – Serum	<b>g/dL</b>
Units – Urine	
Number of needle washes	<b>2/2</b>
Number of cuvette washes	<b>2</b>
Dynamic Blank	<b>Inactive</b>
Reagent Blank	<b>Every Run</b>
Reagent Limit (mABS)	<b>1500</b>
Curve Acceptance (%)	<b>100</b>
Instrument Factor	<b>0.000</b>
Shift	<b>0.000</b>

SERUM	
Name	<b>PROT</b>
Sample (µL)	<b>6</b>
Pré-Dilution	<b>1.00</b>
<u>Dilution</u>	
Factor	<b>2.00</b>
Test Limit	<b>14.0</b>
Initial ABS	<b>@</b>
Final ABS	<b>@</b>
Max ABS Delta	<b>500</b>
<u>Normal Values</u>	
Man	<b>6.0 – 8.0</b>
Woman	<b>6.0 – 8.0</b>
Child	<b>6.0 – 8.0</b>
Calibrator	
Std. 1	<b>0</b>
Std. 2	<b>#</b>

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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### Instrument Settings

Test Methodology	
Kind of Process	<b>End Point</b>
Type	<b>Linear</b>
Filters	<b>510 / 630</b>
Reaction Direction	<b>Increase</b>
Reagent #1 (µL)	<b>300</b>
Reagent #2 (µL)	<b>-</b>
Sample Starter	<b>Inactive</b>
Delay Time (sec.)	<b>0</b>
Incubation Time (sec.)	<b>300</b>
Reading Time (sec.)	<b>20</b>

SERUM	
Name	<b>TRIG</b>
Sample (µL)	<b>3</b>
Pré-Dilution	<b>1.00</b>
<u>Dilution</u>	
Factor	<b>5.00</b>
Test Limit	<b>1100</b>
Initial ABS	<b>@</b>
Final ABS	<b>@</b>
Max ABS Delta	<b>99999</b>
<u>Normal Values</u>	
Man	<b>60 - 150</b>
Woman	<b>60 - 150</b>
Child	<b>60 - 150</b>
Calibrator	
Std. 1	<b>0</b>
Std. 2	<b>#</b>

Units – Serum	<b>mg/dL</b>
Units – Urine	
Number of needle washes	<b>1</b>
Number of cuvette washes	<b>1</b>
Dynamic Blank	<b>Inactive</b>
Reagent Blank	<b>Every Run</b>
Reagent Limit (mABS)	<b>500</b>
Curve Acceptance (%)	<b>20</b>
Instrument Factor	<b>0.000</b>
Shift	<b>0.000</b>

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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## Instrument Settings

Test Methodology	
Kind of Process	<b>End Point</b>
Type	<b>Linear</b>
Filters	<b>578</b>
Reaction Direction	<b>Increase</b>
Reagent #1 (µL)	<b>300</b>
Reagent #2 (µL)	<b>300</b>
Sample Starter	<b>Inactive</b>
Delay Time (sec.)	<b>0</b>
Incubation Time (sec.)	<b>300 / 10</b>
Reading Time (sec.)	<b>300</b>

Units – Serum	<b>mg/dL</b>
Units – Urine	
Number of needle washes	<b>3/3</b>
Number of cuvette washes	<b>3</b>
Dynamic Blank	<b>Inactive</b>
Reagent Blank	<b>Every Run</b>
Reagent Limit (mABS)	<b>500</b>
Curve Acceptance (%)	<b>20</b>
Instrument Factor	<b>0.000</b>
Shift	<b>0.000</b>

SERUM	
Name	<b>UREA UV</b>
Sample (µL)	<b>3</b>
Pré-Dilution	<b>1.00</b>
<u>Dilution</u>	
Factor	<b>1.00</b>
Test Limit	<b>300</b>
Initial ABS	<b>@</b>
Final ABS	<b>@</b>
Max ABS Delta	<b>99999</b>
<u>Normal Values</u>	
Man	<b>15 – 40</b>
Woman	<b>15 – 40</b>
Child	<b>-</b>
Calibrator	
Std. 1	<b>0</b>
Std. 2	<b>#</b>

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

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## Instrument Settings

Test Methodology	
Kind of Process	<b>Fixed time</b>
Type	<b>Linear</b>
Filters	<b>340 / 630</b>
Reaction Direction	<b>Decrease</b>
Reagent #1 (µL)	<b>240</b>
Reagent #2 (µL)	<b>60</b>
Sample Starter	<b>Active</b>
Delay Time (sec.)	<b>0 / 120</b>
Incubation Time (sec.)	<b>60</b>
Reading Time (sec.)	<b>60</b>

Units – Serum	<b>mg/dL</b>
Units – Urine	
Number of needle washes	<b>1/1</b>
Number of cuvette washes	<b>1</b>
Dynamic Blank	<b>Inactive</b>
Reagent Blank	<b>Every Run</b>
Reagent Limit (mABS)	<b>900</b>
Curve Acceptance (%)	<b>20</b>
Instrument Factor	<b>0.000</b>
Shift	<b>0.000</b>

SERUM	
Name	<b>UREA UV</b>
Sample (µL)	<b>3</b>
Pré-Dilution	<b>1.00</b>
<u>Dilution</u>	
Factor	<b>10.00</b>
Test Limit	<b>300</b>
Initial ABS	<b>@</b>
Final ABS	<b>@</b>
Max ABS Delta	<b>1000</b>
<u>Normal Values</u>	
Man	<b>15 – 40</b>
Woman	<b>15 – 40</b>
Child	<b>-</b>
Calibrator	
Std. 1	<b>0</b>
Std. 2	<b>#</b>

URINE	
Name	
Sample (µL)	
Pré-Dilution	
<u>Dilution</u>	
Factor	
Test Limit	
Initial ABS	
Final ABS	
Max ABS Delta	
<u>Normal Values</u>	
Man	
Woman	
Child	

EDIÇÃO: 09/11