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TEST NAME	AU
DECIMAL	1
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	8
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	300
LINEARITY LIMIT	25.0
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	510
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	300
INCUBATION 2	0
LAG PHASE	2
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD *
FACTOR	0.0

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TEST NAME	AEO
DECIMAL	0
MEASURE UNIT	UI/mL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	4
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2.000
LINEARITY LIMIT	800
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	546
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.2
MIX 2	0.00
INCUBATION 1	180
INCUBATION 2	0
LAG PHASE	2
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD *
FACTOR	

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TEST NAME	ALB
DECIMAL	1
MEASURE UNIT	g/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	4
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2.000
LINEARITY LIMIT	7.0
CONTAMINANTING	YES
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	620
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	120
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD *
FACTOR	0

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TEST NAME	ALT
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	40
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	900
MAXIMUM	2000
LINEARITY LIMIT	250
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	40
INCUBATION 2	0
LAG PHASE	20
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	- 1910

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TEST NAME	AMILASE
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	8
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	1300
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	405
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.2
MIX 2	0.00
INCUBATION 1	40
INCUBATION 2	0
LAG PHASE	20
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	3292

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TEST NAME	AMIL
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	
SERUM	4
PLASMA	4
URINE	8
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2400
LINEARITY LIMIT	7400
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	405
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.40
MIX 2	0.00
INCUBATION 1	52
INCUBATION 2	0
LAG PHASE	10
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	24830

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TEST NAME	AST
DECIMAL	1
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	40
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	900
MAXIMUM	2000
LINEARITY LIMIT	250
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	40
INCUBATION 2	0
LAG PHASE	20
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	- 1910

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TEST NAME	BD
DECIMAL	1
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	
SERUM	40
PLASMA	40
URINE	
REAGENT VOLUME (µL)	
REAGENT 1	320
REAGENT 2	80
ABSORBANCE RANGE (mABS)	
MINIMUM	-0.100
MAXIMUM	1000
LINEARITY LIMIT	15
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	YES
FILTERS (nm)	
FILTER 1	546
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.0
MIX 2	0.40
INCUBATION 1	5
INCUBATION 2	300
LAG PHASE	2
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	

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TEST NAME	BT
DECIMAL	1
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	
SERUM	40
PLASMA	40
URINE	
REAGENT VOLUME (µL)	
REAGENT 1	320
REAGENT 2	80
ABSORBANCE RANGE (mABS)	
MINIMUM	-0.100
MAXIMUM	1000
LINEARITY LIMIT	15
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	YES
FILTERS (nm)	
FILTER 1	546
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.0
MIX 2	0.40
INCUBATION 1	5
INCUBATION 2	120
LAG PHASE	2
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	

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TEST NAME	BD
DECIMAL	1
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	
SERUM	24
PLASMA	24
URINE	24
REAGENT VOLUME (µL)	
REAGENT 1	200
REAGENT 2	200
ABSORBANCE RANGE (mABS)	
MINIMUM	-0.100
MAXIMUM	1000
LINEARITY LIMIT	15
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	YES
FILTERS (nm)	
FILTER 1	535
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.0
MIX 2	0.40
INCUBATION 1	5
INCUBATION 2	120
LAG PHASE	2
MEASURE	
MEASUREMENT TYPE	
CALIBRATE	
FACTOR	19 *

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TEST NAME	BT
DECIMAL	1
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	
SERUM	24
PLASMA	24
URINE	24
REAGENT VOLUME (µL)	
REAGENT 1	200
REAGENT 2	200
ABSORBANCE RANGE (mABS)	
MINIMUM	-0.100
MAXIMUM	1000
LINEARITY LIMIT	15
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	YES
FILTERS (nm)	
FILTER 1	535
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.0
MIX 2	0.40
INCUBATION 1	5
INCUBATION 2	120
LAG PHASE	2
MEASURE	
MEASUREMENT TYPE	
CALIBRATE	
FACTOR	19 *

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TEST NAME	CA
DECIMAL	1
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	4
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2.000
LINEARITY LIMIT	15.0
CONTAMINATING	YES
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	620
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	180
INCUBATION 2	0
LAG PHASE	2
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD *
FACTOR	0

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TEST NAME	COLT
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	4
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	300
LINEARITY LIMIT	1000
CONTAMINANTING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	510
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	300
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD *
FACTOR	0

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TEST NAME	HDL
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	40
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	300
LINEARITY LIMIT	150
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	510
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	600
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD *
FACTOR	0

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TEST NAME	LDL
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	8
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	500
LINEARITY LIMIT	500
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	510
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	580
INCUBATION 2	0
LAG PHASE	2
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD *
FACTOR	0

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TEST NAME	CHE
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	6
REAGENT VOLUME (µL)	
REAGENT 1	360
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	25000
CONTAMINATING	YES
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	405
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.0
MIX 2	0.0
INCUBATION 1	90
INCUBATION 2	0
LAG PHASE	10
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	- 65804

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TEST NAME	CK
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	10
REAGENT VOLUME (µL)	
REAGENT 1	500
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	1.000
LINEARITY LIMIT	1.500
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	170
INCUBATION 2	0
LAG PHASE	15
MEASURE	50
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	8284

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TEST NAME	CK MB
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	16
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0.000
MAXIMUM	2.000
LINEARITY LIMIT	600
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	300
INCUBATION 2	0
LAG PHASE	20
MEASURE	80
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	8254

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TEST NAME	CREAT
DECIMAL	2
MEASURE UNIT	mg/dL
REACTION TYPE	FIXED TIME
SAMPLE VOLUME (µL)	40
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2.000
LINEARITY LIMIT	20.0
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	510
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	5
INCUBATION 2	0
LAG PHASE	20
MEASURE	60
MEASUREMENT TYPE	
CALIBRATE	STANDARD *
FACTOR	0

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TEST NAME	CREA
DECIMAL	1
MEASURE UNIT	mg/dL
REACTION TYPE	FIXED TIME
SAMPLE VOLUME (µL)	
SERUM	40
PLASMA	40
URINE (dil. 1/50)	40
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2400
LINEARITY LIMIT	10
CONTAMINATING	YES
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	492
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.40
MIX 2	0.00
INCUBATION 1	60
INCUBATION 2	0
LAG PHASE	10
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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TEST NAME	FERRO
DECIMAL	0
MEASURE UNIT	µg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	80
REAGENT VOLUME (µL)	
REAGENT 1	320
REAGENT 2	80
ABSORBANCE RANGE (mABS)	
MINIMUM	0.0
MAXIMUM	2.000
LINEARITY LIMIT	625
CONTAMINATING	YES
REAGENT BLANKING	YES
DIFFERENTIAL	YES
FILTERS (nm)	
FILTER 1	546
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.0
MIX 2	0.2
INCUBATION 1	5
INCUBATION 2	300
LAG PHASE	2
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD *
NORMAL RANGE	55-155

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TEST NAME	FER
DECIMAL	0
MEASURE UNIT	µg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	
SERUM	80
PLASMA	80
URINE	80
REAGENT VOLUME (µL)	
REAGENT 1	200
REAGENT 2	200
ABSORBANCE RANGE (mABS)	
MINIMUM	-0.100
MAXIMUM	1000
LINEARITY LIMIT	600
CONTAMINATING	YES
REAGENT BLANKING	YES
DIFFERENTIAL	YES
FILTERS (nm)	
FILTER 1	546
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0
MIX 2	0.4
INCUBATION 1	5
INCUBATION 2	300
LAG PHASE	2
MEASURE	
MEASUREMENT TYPE	
CALIBRATE	
FACTOR	1380 *

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TEST NAME	FAC
DECIMAL	1
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	40
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	1200
LINEARITY LIMIT	150
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	405
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	240
INCUBATION 2	0
LAG PHASE	10
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	844

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TEST NAME	FALC
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	8
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	1200
LINEARITY LIMIT	1200
CONTAMINATING	YES
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	405
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	60
INCUBATION 2	0
LAG PHASE	20
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	2764

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TEST NAME	FOSF
DECIMAL	1
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	
SERUM	4
PLASMA	4
URINE	-
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2400
LINEARITY LIMIT	20
CONTAMINATING	YES
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.40
MIX 2	0.00
INCUBATION 1	300
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDAR D
FACTOR	0

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TEST NAME	FRUT
DECIMAL	2
MEASURE UNIT	mmol/L
REACTION TYPE	FIXED TIME
SAMPLE VOLUME (µL)	20
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	500
LINEARITY LIMIT	7.0
CONTAMINATING	YES
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	535
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	600
INCUBATION 2	0
LAG PHASE	10
MEASURE	300
MEASUREMENT TYPE	
CALIBRATE	STANDARD *
FACTOR	0

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TEST NAME	GGT
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	30
REAGENT VOLUME (µL)	
REAGENT 1	300
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2.000
LINEARITY LIMIT	300
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	405
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	60
INCUBATION 2	0
LAG PHASE	20
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	1160

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TEST NAME	GGT
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	
SERUM	40
PLASMA	40
URINE	-
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	300
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	405
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.60
MIX 2	0.00
INCUBATION 1	52
INCUBATION 2	0
LAG PHASE	10
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	1158

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TEST NAME	GLICOSE
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	4
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	300
LINEARITY LIMIT	500
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	510
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	250
INCUBATION 2	0
LAG PHASE	2
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD *
FACTOR	0

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TEST NAME	HDL DIRETO
DECIMAL	1
MEASURE UNIT	mg/dL
REACTION TYPE	FIXED TIME
SAMPLE VOLUME (µL)	3
REAGENT VOLUME (µL)	
REAGENT 1	300
REAGENT 2	100
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	200
CONTAMINANTING	YES
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	546
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.20
INCUBATION 1	300
INCUBATION 2	5
LAG PHASE	3
MEASURE	120
MEASUREMENT TYPE	
CALIBRATE	STANDARD *
FACTOR	

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TEST NAME	LDL DIRETO
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	3
REAGENT VOLUME (µL)	
REAGENT 1	375
REAGENT 2	125
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	450
CONTAMINANTING	YES
REAGENT BLANKING	YES
DIFFERENTIAL	YES
FILTERS (nm)	
FILTER 1	546
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.20
INCUBATION 1	300
INCUBATION 2	300
LAG PHASE	2
MEASURE	2
MEASUREMENT TYPE	
CALIBRATE	STANDARD *
FACTOR	

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TEST NAME	LDH
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	8
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	800
MAXIMUM	2000
LINEARITY LIMIT	1250
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	40
INCUBATION 2	0
LAG PHASE	20
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	- 8095

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TEST NAME	MG
DECIMAL	2
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	4
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2400
LINEARITY LIMIT	4.0
CONTAMINANTING	YES
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	535
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	180
INCUBATION 2	0
LAG PHASE	2
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD *
FACTOR	0

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TEST NAME	ALB
DECIMAL	1
MEASURE UNIT	mg/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	3
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	600
LINEARITY LIMIT	130
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	546
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	5
INCUBATION 2	0
LAG PHASE	90
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	STANDARD *
FACTOR	0.0

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TEST NAME	PT
DECIMAL	1
MEASURE UNIT	g/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	8
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	300
LINEARITY LIMIT	20
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	535
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	240
INCUBATION 2	0
LAG PHASE	2
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD *
FACTOR	0

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TEST NAME	PROTUR
DECIMAL	1
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	20
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	160
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	620
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	550
INCUBATION 2	0
LAG PHASE	2
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD *
FACTOR	0

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TEST NAME	TRIG
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	4
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	300
LINEARITY LIMIT	600
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	510
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	300
INCUBATION 2	0
LAG PHASE	2
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD *
FACTOR	0

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TEST NAME	UREUV
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	FIXED TIME
SAMPLE VOLUME (µL)	3
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	1000
MAXIMUM	2000
LINEARITY LIMIT	300
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	20
INCUBATION 2	0
LAG PHASE	10
MEASURE	60
MEASUREMENT TYPE	
CALIBRATE	STANDARD *
FACTOR	0

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TEST NAME	UREIA
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	3
REAGENT VOLUME (µL)	
REAGENT 1	300
REAGENT 2	300
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	300
CONTAMINATING	YES
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	620
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.20
INCUBATION 1	300
INCUBATION 2	300
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD *
FACTOR	

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© AIRONE 200 é marca registrada de seus proprietários.